


GAMMA

COSY Examples



Author: Scott A. Smith
Date: August 19, 1999

Table of Contents

1	<i>COSY</i>	6
1.1	Introduction	6
1.2	COSY- NMRi, Ideal Pulses	7
1.2.1	Description	7
1.2.2	Program	8
1.2.3	Discussion	8
1.2.4	Spin System	9
1.2.5	Results	9
1.3	COSY- Felix, Ideal Pulses	11
1.3.1	Description	11
1.3.2	Program	12
1.3.3	Discussion	12
1.3.4	Spin System	13
1.3.5	Results	13
1.4	COSY with TPPI	16
1.4.1	Description	16
1.4.2	Program	17
1.4.3	Discussion	17
1.4.4	Spin System	18
1.4.5	Results	18
1.5	COSY with RuSH Method	21
1.5.1	Description	21
1.5.2	Program	22
1.5.3	Discussion	23
1.5.4	Felix Processing	24
1.5.5	Spin System	28
1.5.6	Workup	29
1.5.7	Results	32
1.6	COSY with real pulses	33
1.6.1	Description	33
1.6.2	Program	34
1.6.3	Results	35
2	<i>E-COSY</i>	36
2.1	Introduction	36
2.2	E-COSY Examples	38
2.3	E-COSY with Superoperators	39
2.3.1	Description	39
2.3.2	Program	45
2.3.3	Discussion	46
2.3.4	Spin System	47

2.3.5	Workup	47
2.3.6	Results	49
2.3.7	Spin System	50
2.3.8	Workup	51
2.3.9	Results	52
2.4	Alternative E-COSY	53
2.4.1	Description	53
2.4.2	Program	56
2.4.3	Discussion	57
2.4.4	Spin System	58
2.4.5	Workup	58
2.5	Complimentary E-COSY with Superoperators	59
2.5.1	Description	59
2.5.2	Program	60
2.5.3	Discussion	62
2.5.4	Spin System	63
2.5.5	Workup	63
2.5.6	Results	64
2.6	Heteronuclear E-COSY	65
2.6.1	Description	65
2.6.2	Program	67
2.6.3	Discussion	68
2.6.4	Spin Systems	70
2.6.5	Results	71
2.7	E-COSY Phase Cycling	72
2.7.1	Description	72
2.7.2	Program	73
2.7.3	Discussion	74
2.7.4	Results	76
2.8	E-COSY with Relaxation	77
2.8.1	Description	77
2.8.2	Program	80
2.8.3	Discussion	81
2.8.4	Spin System	83
2.8.5	Workup	84
2.8.6	Results	86
2.9	E-COSY with Phase Cycling	87
2.9.1	Description	87
3	<i>MQF-COSY</i>	92
3.1	Introduction	92
3.2	MQF-COSY with Superoperators & Felix Output	93
3.2.1	Description	93
3.2.2	Program	95

3.2.2	Discussion	96
3.2.3	Example Spin System	97
3.2.4	Workup	98
3.2.5	Results	100
4	<i>NMR Mathematical Description</i>	<i>101</i>
4.1	Overview	101
4.2	Liouville Equation in the Laboratory Frame	101
4.3	Liouville Equation, Interaction Representation	101
4.4	Liouville Equation in the Rotating Frame	102
4.5	Liouville Equation, Multiple Rotating Frame	104
4.6	Liouville Equation Solution, Static Hamiltonian	105
4.7	Example Source Codes	107

1 COSY

1.1 Introduction

The pulse sequence below is the simplest of the COSY experiments.




Figure 0-1 A simple COSY pulse sequence .

The first pulse generates magnetization in the xy-plane. Subsequently, the magnetization is frequency labeled during t_1 and the final pulse produces the observable signal. The chemical shifts of the spin system show up as frequencies during the acquisition and these are modulated by the frequency labeling during t_1 to all spins with which they are spin coupled.

The following set of programs which simulate COSY-experiments are discussed in this document:

Table 1: COSY Example Programs

Example	Page	Pulse ^a	Relaxation	System	Workup ^b
1. COSYNMRi	36	Ideal	No	Homonuclear	NMR2
2. COSYFelix	46	Ideal	No	Homonuclear	Felix
3. COSYTPPI	54	Ideal	No	Homonuclear	Felix/TPPI
4. COSYRuSH	61	Real	No	Heteronuclear	Felix

a. "Ideal" pulses are "infinitely" short so no relaxation effects can be considered during such sequence steps.

b. Workups labeled "Open" allow the user to choose the output type while the program is running. The output types currently available in GAMMA are FrameMaker, Felix, NMRi, and MATLAB

Please note that these examples are meant as simple tutorials in getting COSY programs built. They are typically NOT the most elegant nor the most versatile. However they should provide a wide basis on which robust simulation programs can be built. To obtain these and more COSY programs see the GAMMA WWW site at <http://gamma.magnet.fsu.edu>.

1.2 COSY- NMRI, Ideal Pulses

1.2.1 Description

This example performs a simple COSY simulation. The pulse sequence diagram is repeated with density matrix labels used to coincide with variable names in the program code.




Figure 0-2 A simple COSY pulse sequence. Density operators label key points in simulation of this experiment.

The code itself will not be the most efficient in producing a COSY simulation, efficiency was sacrificed (slightly) for code clarity. No relaxation effects will be considered and the pulses taken as ideal. The resulting output file is produced in NMRI format and subsequently processed with that program. The size of the working data set is set for 1K x 1K and the spin system only has three protons. Although the run time is in minutes, although the data workup is quite lengthy.

1.2.2 Program

```

/* cosy1.cc *****-c++-*/
**
**           Example program for the GAMMA Library
**
** This program simulates a basic COSY experiment.
** No phase cycling is performed and all input are
** kept positive to avoid the need for quadrature
** detection in t1.
**
***** */

#include <gamma.h>
main()
//           DEFINE SYSTEM & NMR PARAMETERS

const double dt1 = 0.001;           // t1 time increment
const double dt2 = 0.001;           // t2 time increment
const int t1pts = 1024;              // points on t1 axis
const int t2pts = 1024;              // points on t2 axis
spin_system sys;                    // define the system, read in
sys.read("cosy.1.sys");              // from disk

//           SET UP NECESSARY VARIABLES

block_1D tmp(t2pts);                // 1D-data block storage
block_2D data(t1pts,t2pts);          // 2D-data matrix storage
gen_op H = Hcs(sys)+ HJw(sys);       // Hamiltonian, weak coupling
gen_op detect = Fm(sys);             // F- for detection operator
gen_op sigma0, sigma1, sigma;        // working density matrices

//           APPLY PULSE SEQUENCE

sigma0 = sigma_eq(sys);              // equilibrium density matrix
sigma1 = lpuls(sys, sigma0, 90);      // apply first 90 y-pulse
for (int t1=0; t1<t1pts; t1++)        // loop over t1 points
{
    sigma=evolve(sigma1, H, t1*dt1);  // evolution during t1
    sigma=lpuls(sys, sigma, 90);       // apply second 90 y-pulse
    FID(sigma,detect,H,dt2,t2pts,tmp); // acquisition
    data.put_block(t1, 0, tmp);        // store it in a 2D block
}
NMRi("cosy.1.dat", data);            // output for NMRi workup
}

```

1.2.3 Discussion

Define System and NMR Parameters: The first four lines specify the dwell times and number of points to use on the t_1 and t_2 axes. The spin system used was a three spin homonuclear system with weak coupling, read in from disk file cosy1.sys (see next section).

Set up the Necessary Variables: For this simulation two data blocks are allocated. The first is a 1D-block called tmp for temporary FID storage. The second is a 2D-block called data into which each 1D block is successively placed. The Hamiltonian is set for weak scalar coupling (under the secular approximation). The detector is set to $F_- = F_x - iF_y$ so as to produce a complex data set along the t_2 axis. Three density matrices are specified, coinciding with the three matrices in the previous pulse sequence figure.

Apply the Pulse Sequence: The program sets the initial density matrix, σ_0 , to the equilibrium density matrix. A ninety degree pulse is then applied to it to produce σ_1 . As can be seen from the applied pulse sequence, the density matrix σ_1 is constant for all mixing times. It's computation is thus performed prior to beginning any looping over t_1 times. Following this the looping begins over all t_1 points. The first loop step is to evolve σ_1 for the t_1 time to produce $\sigma(t_1)$. The second 90 pulse of the sequence is then applied and the xy-magnetization determined. This is an acquisition step which fills the data block tmp up with t2pts FID points. The final step is to take that individual 1D-block and insert it into a 2D-block, data. Then data is the 2D-data set, output on the last step of the program in NMRi format.

1.2.4 Spin System

The spin system file, cosy.1.sys, is reproduced below.

SysName	(2) : COSY3	- Name of the Spin System
NSpins	(0) : 3	- Number of Spins in the System
Iso(0)	(2) : 1H	- Spin Isotope Type
Iso(1)	(2) : 1H	- Spin Isotope Type
Iso(2)	(2) : 1H	- Spin Isotope Type
v(0)	(1) : 200.0	- Chemical Shifts in Hz
v(1)	(1) : 100.0	- Chemical Shifts in Hz
v(2)	(1) : 40.0	- Chemical Shifts in Hz
J(0,1)	(1) : 10.0	- Coupling Constants in Hz
J(1,2)	(1) : 15.0	- Coupling Constants in Hz
Omega	(1) : 400	- Spectrometer Frequency in MHz (1H based)

1.2.5 Results

The 2D-data file produced was processed with NMRi, specifically NMR2^{1,2}. The data set was manipulated within NMR2 by a link (macro resulting from the linking) of the following commands:

SP	apodize in t2 with Sine bell Phased function (0.5, 1, 6)
FT	Fourier Transformation in t2
TP	Transpose the entire data set
EI	set Imaginary part to 0 since it is redundant with real data
SP	apodize in t1 with Sine bell Phased function (0.5, 1, 6)
FT	Fourier transformation in t1
TP	Transpose the entire data set back so F2 is horizontal

Note that the comments to the right are only for description of the two letter commands used to create the link. The link is produced in the following manner.

1. CM - go into Command Mode
 - 1a.) M- Modify/create link.

1a.-1)	- Supply a link name.
1a.-2)	- Type two letter link commands (above).
1a.-3)	- Hit return with no command to end the link.
 - 1b.) X- Execute (in a dry run first) the link

1b-1)	- SP parameters are 0.5, 1, 6
1b-2)	- EI parameter is 0
1b-3)	- SP parameters are 0.5, 1, 6 (same as t2 axis)
2. DM - go into Display Mode. (Note all DM commands not shown in window)
 - 1a.) P- Define contour levels.

1. NMR1 and NMR2 (jointly called NMRi) are products of New Methods Research, Inc. in Syracuse, New York. Phone: (315) 424-0329. These programs are designed specifically for the workup of NMR data.
2. This example produces a 1K by 1K complex output file. Data processing will be time consuming and limited by computer system memory and/or disk space availability. Switching from 1024 to 256 in the program, both time and memory requirements for the simulation and data processing will be greatly reduced.

The COSY spectrum from simulation cosy.1. This is the NMR2 output format

```
Parameters PLP.H01 - NMR2 3.95
Data File: /data4/sosi/shit.trf

1st Dim Full SW: 1000 Hz
2nd Dim Full SW: 1000 Hz
F2 Obs Freq: 400.0 MHz
F1 Obs Freq: 400.0 Mhz
Lock Freq: 0.0 Mhz
Acquisitions: 0

Apodization Function SP.
->F2 Q1 = 0.5
->F2 Q2 = 1.0
->F2 Q3 = 6.0

F2 Zero Fills: 0
F2 Phase: (0.00,0.00)

Size: 1024
Spectra: 1024
Data Type: quadrature
Data are not transposed.
Mode: Phase Sens. TPPI.
FDATA(512): 0.0

Apodization Function SP.
->F1 Q1 = 0.5
->F1 Q2 = 1.0
->F1 Q3 = 6.0

F1 Zero Fills: 0
F1 Phase: (0.00,0.00)
```




Figure 0-3 Simulated COSY spectrum output from NMR2, the result of COSY1.cc. The plot was placed into eps format via HPGL output from NMR2 & the filter hpgltoeps provided by FrameMaker (used for this document)

1.3 COSY- Felix, Ideal Pulses

1.3.1 Description

This example again performs a basic COSY simulation. The pulse sequence diagram is repeated with density matrix labels used to coincide with variable names in the program code.




Figure 0-4 A simple COSY pulse sequence. Density operators label key points in simulation of this experiment.

The code itself is written to be a little more efficient than in the previous example. No relaxation effects will be considered and the pulses taken as ideal. The resulting output file is produced in Felix¹ format and subsequently processed with that program. There are two differences of note between this example and the previous one. First the data processing here provides a clear example of how to work up simulated data with the Felix program (rather than NMR2 used previously). Secondly, use is made of a pulse propagator rather than recomputing the pulses for each step and each loop of the pulse sequence. Mathematically this is given by

$$\sigma_1 = U \sigma_{eq} U^{-1} \quad \text{and} \quad \sigma'(t_1) = U \sigma(t_1) U^{-1}$$

where U is the 90-y pulse propagator given by

$$U = \exp[-(iF_y \pi)/2] .$$

1. Felix is an NMR data processing program by Hare Research, Inc.

1.3.2 Program

```

/* cosy2.cc *****-c++- */
**
**           Example program for the GAMMA Library
**
** This program simulates a basic COSY experiment.
** No phase cycling is performed and all input are
** kept positive to avoid the need for quadrature
** detection in t1.
**
***** */

#include <gamma.h>
main()
{
//           DEFINE SYSTEM & NMR PARAMETERS
const double dt1 = 0.0015;           // t1 time increment
const double dt2 = 0.0015;           // t2 time increment
const int t1pts = 512;                // points on t1 axis
const int t2pts = 512;                // points on t2 axis
spin_system sys;                     // define the system, read in
sys.read("cosy.2.sys");               // from disk

//           SET UP NECESSARY VARIABLES
block_1D data(t2pts);                // 1D-data block storage
gen_op H = Hcs(sys)+ HJw(sys);        // Hamiltonian, weak coupling
gen_op detect = Fm(sys);              // F- for detection operator
gen_op y_pulse = lypuls_U(sys, 90);    // 90 y-pulse propagator
gen_op sigma0, sigma1, sigma;         // working density matrices

//           APPLY PULSE SEQUENCE
sigma0 = sigma_eq(sys);               // equilibrium density matrix
sigma1 = evolve(sigma0, y_pulse);      // apply first 90 y-pulse
File fp;
fp.open("felix.dat", io_writeonly, a_create);
for (int t1=0; t1<t1pts; t1++)        // loop over t1 points
{
    sigma = evolve(sigma1, H, t1*dt1); // evolution during t1
    evolve_ip(sigma, y_pulse);          // apply second 90 y-pulse
    FID(sigma, detect, H, dt2, t2pts, data); // acquisition
    Felix(fp, data);
}
fp.close();
}

```

1.3.3 Discussion

Define System and NMR Parameters: The first four lines specify the dwell times and number of points to use on the t_1 and t_2 axes. The spin system used was a three spin homonuclear system with weak coupling, read in from disk file cosy.2.sys (see next section).

Set up the Necessary Variables: For this simulation only one data block is allocated (unlike the two in the previous example). The Hamiltonian is set for weak scalar coupling (under the secular approximation). The detector is set to $F_- = F_x - iF_y$ so as to produce a complex data set along the t_2 axis. The next operator is a propagator for a y-pulse of 90 degrees. Since it is repeatedly used in the pulse sequence it is computationally efficient to produce it outside of the loop. Three density matrices are specified, coinciding with the three matrices in the previous pulse sequence figure.

Apply the Pulse Sequence: The program sets the initial density matrix, σ_0 , to the equilibrium density matrix. A ninety degree pulse is applied using the y-pulse propagator to produce σ_1 . As can be seen from the applied pulse sequence, the density matrix σ_1 is constant for all mixing times. It's computation is thus performed prior to beginning any looping over t_1 times. Following this the looping begins over all t_1 points. Just before the loop a file is opened called "felix.dat" for the spectral output. The first loop step is to evolve σ_1 for the t_1 time to produce $\sigma(t_1)$. The second 90 pulse of the sequence is then applied and the xy-magnetization determined. This is an acquisition step which fills the data block data up with t2pts FID points. The final step in the loop is to write this block out to the file "felix.dat" in Felix format. The last line of the program simply closes the file.

1.3.4 Spin System

The spin system file, cosy2.sys, is virtually identical to the one used in the previous simulation. It is reproduced below.

SysName	(2) : COSY3	- Name of the Spin System
NSpins	(0) : 3	- Number of Spins in the System
Iso(0)	(2) : 1H	- Spin Isotope Type
Iso(1)	(2) : 1H	- Spin Isotope Type
Iso(2)	(2) : 1H	- Spin Isotope Type
v(0)	(1) : 200.0	- Chemical Shifts in Hz
v(1)	(1) : 100.0	- Chemical Shifts in Hz
v(2)	(1) : 40.0	- Chemical Shifts in Hz
J(0,1)	(1) : 10.0	- Coupling Constants in Hz
J(0,2)	(1) : 0.0	- Coupling Constants in Hz
J(1,2)	(1) : 15.0	- Coupling Constants in Hz
Omega	(1) : 400	- Spectrometer Frequency in MHz (1H based)

1.3.5 Results

Workup of the simulated data was performed in Felix. A text editor was used to create the following Felix macro file, called cosy.2.mac (the comments off to the side are not part of the macro) -

cl	! insure any data files are closed
cmx	! insure any matrix files are closed
get 'Enter matrix name: ' mname	! ask user for matrix
mat &mname write	! open matrix with write ability
get 'Enter FID name: ' dname	! ask user for fid file name
lb 2	! set line broadening
for row 1 512	! loop through all 512 rows
re &dname	! read FID from file
em	! apodize
zf 1024	! zero fill
ft	! fourier transform
rev	! reverse spectrum so frequencies proper
si 512	! reduce the size back to 512
sto 0 &row	! store block in matrix row
ty row = &row	! type row is complete
next	! get next row
for col 1 512	! begin column processing
loa &col 0	! retrieve column from matrix
em	! apodize
zi	! zero imaginaries (real data on this axis)
zf 1024	! zero fill
red	! reduced from complex to real
rft	! real fourier transform
rev	! reverse so correct frequencies
sto &col 0	! store column back into matrix
ty col = &col	! type that column is complete
next	! get next column
end	

Once Felix has been started, the first task is to build a matrix into which the data is to be placed.

The command (in Felix) for this is¹

```
bld cosy 2 512 512 1
```

This creates a matrix files called cosy.mat which is 512x512 complex. Now the macro listed previously can be run which will read the data in and process it. The command for this (still in Felix) is

```
ex cosy.2.mac
```

where cosy.2.mac is the name of the file containing the macro². The macro will ask you to input the matrix file name (cosy as named by the bld command) and the FID file name (felix as determined by the program when creating the output). Once this macro has been executed in Felix, the COSY spectrum is ready to be displayed. Commands to produce the contour plot on the screen are³ as follows.

lvl 8e-4	! set contour level low (~10 ⁸ lower than real data)
cpn 1	! both positive and negative contours
nl 5	! five contour levels
cyc 2	! even pen cycle to show positive and negatives
clm 2.5	! contours increment 250 percent
cp	! draw the contour plot to the screen

To get the axes to have the proper scaling one needs to use the rmx command⁴. The spectral width on both axes should be 333.33 Hz. To generate the hpgl contour plot file, the Felix commands are (once cp produces the screen plot)

hdv felix.hpgl	! hardcopy device is file felix.hpgl
hpm 32	! hardcopy plot mode is hpgl
hcp	! produce the plot

This output file can be sent to a plotter which understands HPGL. As a last step for this example the felix output file, felix.hpgl, was converted to MIF format⁵ in order to import the spectrum into

-
1. The matrix used is 512 x 512 complex. Currently the Felix documentation is contradictory in its statements on how to accomplish this. Experience says that either a 1 or a 2 at the end of the build statement produces a complex array, only a 0 produces a real array.
 2. Felix may have trouble finding this macro file unless the "pre" command is used to specify which directory the macro files are located in, see the Felix documentation. In Felix, make sure that directory names end with a /. Keep in mind that Felix has problems (at least in UNIX) with capitol letters in filenames, best not to use and upper-case letters in your Felix related filenames. Finally, Felix has trouble with long path names. If your files are buried deep in some subdirectory use a link to shorten the path length down (see the UNIX man pages: man ln).
 3. Early versions 1.0 of Felix seem to prefer that one draw 1D plots before 2D plots or it confuses itself on the plot limits. Apparently loading a row (loa 0 256) and drawing it (dr) then loading a column (loa 130 0) and drawing it prior to the contour plot does something to help Felix figure itself out. This seems not be a problem in later releases.
 4. The command rmx is not fully documented in the Felix 1.0 manual. However there are scattered references to it throughout the examples. The command will query the user for needed parameters.
 5. MIF stands for Maker Interchange Format used by FrameMaker. This conversion was done with a private filter function called Felix2mif and should reside (after GAMMA installation) in the GAMMA directory which contains the cosy example files. A similar program is supposed to be supplied by the manufacturers of FrameMaker in the near future. The program hpgltoeps that is supplied will perform a similar function and enables to spectrum to be imported into documents in encapsulated postscript format.

this FrameMaker document. The spectrum is nearly identical to that in the previous simulation.

3 Spin COSY Simulation, Felix Workup




Figure 0-5 Simulated COSY spectrum output from Felix, the result of COSY2.cc

1.4 COSY with TPPI

1.4.1 Description

This example again performs a simple COSY simulation. In this instance, a TPPI¹ scheme is implemented in order to achieve quadrature detection along the t_1 axis. This is done by phase incrementing the initial pulse by 90 degrees on successive points along t_1 . For clarity, the pulse sequence diagram is repeated with density matrices used to coincide with variable names in the program code.




Figure 0-6 A simple COSY pulse sequence. Density operators label key points in simulation of this experiment.

In this example, rather than incrementing the receiver phase (normal TPPI), for convenience the first pulse phase is incremented. No relaxation effects will be considered and the pulses taken as ideal. The resulting output file is produced in Felix format and subsequently processed with that program.

1. The TPPI detection scheme is discussed in Derome, *Modern NMR Techniques for Chemistry Research*, page 83.

1.4.2 Program

```

/* cosy3.cc *****-c++-*/
**
**          Example program for the GAMMA Library
**
** This program simulates a basic COSY experiment. Phase cycling is
** performed to achieve quadrature detection in t1.
**
***** */
#include <gamma.h>
main ()
{
//          DEFINE SYSTEM & NMR PARAMETERS
const int t1pts = 512;          // points on t1 axis
const int t2pts = 512;          // points on t2 axis
const double dt1 = 0.002;       // t1 increment (SW 500 Hz, TPPI)
const double dt2 = 0.004;       // t2 increment (SW 250 Hz)
spin_system sys;                // define the system, read in
sys.read("cosy.3.sys");         // from disk
//          SET UP NECESSARY VARIABLES
gen_op sigma1[4];               // density matrices, for TPPI
gen_op sigma;                   // working density matrix
gen_op Upulse = lxpuls_U(sys, 90.0); // 90 x-pulse propagator
gen_op detect = Fm(sys);        // F- for detection (0 phase)
gen_op H = Hcs(sys) + HJw(sys); // Hamiltonian, weak coupling
gen_op Udelay = le(sys,0);       // t1 propagator (start = 1)
gen_op Udelay1 = prop(H, dt1);   // t1 propagator increment
block_1D data(t2pts);           // Storage for the FID
//          APPLY PULSE SEQUENCE
gen_op sigma0 = sigma_eq(sys);   // equilib. density mx
sigma1[0] = lxpuls(sys,sigma0, 90.0, 0.0); // Apply a (PI/2)x pulse
sigma1[1] = lxpuls(sys,sigma0, 90.0, 90.0); // Apply a (PI/2)y pulse
sigma1[2] = lxpuls(sys,sigma0, 90.0, 180.0); // Apply a (PI/2)-x pulse
sigma1[3] = lxpuls(sys,sigma0, 90.0, 270.0); // Apply a (PI/2)-y pulse
File fp;
fp.open("felix.dat", io_writeonly, a_create); // Create file,open
for(int t1=0; t1<t1pts; t1++)
{
    sigma = evolve(sigma1[t1%4], Udelay); // TPPI Cycle, 1st pulse
    evolve_ip(sigma, Upulse);             // apply 2nd 90 x-pulse
    FID(sigma, detect, H, dt2, t2pts, data); // acquisition
    Felix(fp, data);                     // write block
    Udelay *= Udelay1;                   // increment t1
}
fp.close();
}

```

1.4.3 Discussion

Define System and NMR Parameters: The first four lines specify the dwell times and number of points to use on the t_1 & t_2 axes. Note that the spectral width is set two times that actually desired along t_1 due to use of TPPI. The spin system used was a slight variation of the three spin homonuclear system used in the first two examples. It is read in from disk file cosy3.sys (see next section).

Set up the Necessary Variables: Initially four density matrices are set up for the four phases to be used in TPPI mode following the first 90 pulse. Another working density matrix is also allocated. Next, a pulse propagator is created for the 90x pulse and will be used inside the loop over t_1 increments. The detector is set to $F_- = F_x - iF_y$ so as to produce a complex data set along the t_2 axis. The Hamiltonian is set for weak scalar coupling (under the secular approximation). Two more propagators are now set up. The first is a working propagator which will propagate the density matrix through the total t_1 time. This is initialized to the identity matrix. The 2nd propagator is that for the evolution over a single t_1 increment. Following this, a data block is allocated for the acquisitions.

Apply the Pulse Sequence: The program sets the initial density matrix, σ_0 , to the equilibrium density matrix. A ninety degree pulse is now applied, four times at different phases to produce four σ_1 matrices corresponding to the four phases. These are calculated outside the loop over t_1 times. Just before the loop a file is opened called "felix.dat" for the spectral output. The looping begins over all t_1 points. The first loop step is to evolve the appropriate σ_1 for the t_1 time to produce $\sigma(t_1, \phi)$. In this instance, modulus arithmetic is utilized to convert the t_1 increment into an integer 0,1,2, or 3 corresponding to the phases 0, 90, 180, and 270. The second 90 pulse of the sequence is then applied and the xy-magnetization determined. This is an acquisition step which fills the data block data up

with t2pts FID points. The next step in the loop is to write this block out to the file “felix.dat” in Felix format. At the end of the loop, the propagator for the t_1 delay is incremented to the total delay time in the loop. The last line of the program simply closes the output file.

1.4.4 Spin System

The spin system file, cosy3.sys, is below. It is frequency shifted from the one used previously.

SysName	(2) : COSY_TPPI	- Name of the Spin System
NSpins	(0) : 3	- Number of Spins in the System
Iso(0)	(2) : 1H	- Spin Isotope Type
Iso(1)	(2) : 1H	- Spin Isotope Type
Iso(2)	(2) : 1H	- Spin Isotope Type
v(0)	(1) : 200.0	- Chemical Shifts in Hz
v(1)	(1) : 100.0	- Chemical Shifts in Hz
v(2)	(1) : 40.0	- Chemical Shifts in Hz
J(0,1)	(1) : 10.0	- Coupling Constants in Hz
J(0,2)	(1) : 0.0	- Coupling Constants in Hz
J(1,2)	(1) : 15.0	- Coupling Constants in Hz
Omega	(1) : 400	- Spectrometer Frequency in MHz (1H based)

1.4.5 Results

Workup of the simulated data was performed in Felix. A text editor was used to create the following Felix macro file, called cosy3.mac (the comments off to the side are not part of the macro) -

cl	! insure any data files are closed
cmx	! insure any matrix files are closed
get 'Enter matrix name: ' mname	! ask user for matrix
mat &mname write	! open matrix with write ability
get 'Enter FID name: ' dname	! ask user for fid file name
lb 4	! set line broadening
ph0 90	! set zero order phase correction
for row 1 512	! loop through all 512 rows
re &dname	! read FID from file
em	! apodize
ft	! fourier transform
ph	! apply phase correction
rev	! reverse spectrum so frequencies proper
si 512	! reduce the size back to 512
sto 0 &row	! store block in matrix row
ty row = &row	! type row is complete
next	! get next row
for col 1 512	! begin column processing
loa &col 0	! retrieve column from matrix
em	! apodize
zf 1024	! zero fill
rft	! real fourier transform
red	! reduced from complex to real
rev	! reverse so correct frequencies
sto &col 0	! store column back into matrix
ty col = &col	! type that column is complete

```
next          ! get next column
end
```

Once Felix has been started, the first task is to build a matrix into which the data is to be placed. The command for this is (in Felix)¹

```
bld cosy 2 512 512 0
```

This creates a matrix files called cosy.mat which is 512x512 real. Now the macro listed previously can be run which will read the data in and process it. The command for this (still in Felix) is

```
ex cosy3
```

where cosy3.mac is the name of the file containing the macro². Once this macro has been executed in Felix, the COSY spectrum is ready to be displayed. Commands to produce the contour plot on the screen are³ -

```
lvl 1.3e-4      ! set contour level very low (10**8 lower than exptl.)
cpn 1           ! both positive and negative contours
nl 5            ! five contour levels
cyc 2           ! even pen cycle to show positive and negatives
cp             ! draw the contour plot to the screen
```

To get the axes to have the proper scaling one needs to use the rmx command⁴. The spectral widths on both axes should be 250 Hz and the 256th point set to 120 Hz. To generate the hpgl contour plot file, the Felix commands are (once cp produces the screen plot)

```
hdv felix.hpgl  ! hardcopy device is file felix.hpgl
hpm 32          ! hardcopy plot mode is hpgl
hcp            ! produce the plot
```

As a last step, the felix output file, felix.hpgl was converted to mif format in order to import the spectrum into this document. The spectrum is nearly identical to that in the previous simulation.

-
1. The matrix used is 512 x 512 real. Currently the Felix documentation is contradictory in its statements on how to accomplish this. Experience says that a 0 at the end of the build statement produces a real array, a 1 or 2 produces a complex array.
 2. Felix may have trouble finding this macro file unless the "pre" command is used to specify which directory the macro files are located in. Keep in mind that Felix has problems (at least in UNIX) with capitol letters in filenames and that directory names need to end with a /.
 3. Early versions 1.0 of Felix seem to prefer that one draw 1D plots before 2D plots or it confuses itself on the plot limits. Apparently loading a row (loa 0 256) and drawing it (dr) then loading a column (loa 130 0) and drawing it prior to the contour plot does something to help Felix figure itself out. This seems not be a problem in later releases.
 4. The command rmx is not fully documented in the Felix 1.0 manual. However there are scattered references to it throughout the examples. The command will query the user for needed parameters.

3 Spin COSY Simulation, TPPI, Felix Workup




Figure 0-7 Simulated COSY spectrum from ccosy3.cc. The system used in this simulation had all frequencies shifted by 120 Hz, effectively placing the carrier in the middle of the spectrum. This value was reinserted in the workup so that the spectrum has the same frequencies as in previous examples.

1.5 COSY with RuSH Method

1.5.1 Description

For our final example of a simple COSY simulation we use the RuSH method¹ in order to achieve quadrature detection along the t_1 axis. To accomplish this we essentially perform two simultaneous simulations (experiments), the second being a phase shifted version of the first. Data from the first experiment is then considered during t_1 processing as the “real” data and that of the second experiment as the “imaginary” data. The pulse sequence diagram to be utilized is shown below with density matrices labeled to coincide with the variable names used in the program code.





Figure 0-8 Two COSY pulse sequences for RUSH. Density operators label key points in simulation.

In this simulation two data sets (labeled S1 and S2 in the figure) are generated. Each contain different phase information and they need to be combined for the quadrature detection in t_1 . They are recombined during spectral workup in the following fashion.

COSY Data Workup for RuSH Quadrature Detection



1. The RuSH method is discussed in Derome, *Modern NMR Techniques for Chemistry Research*, page 201.

1.5.2 Program

```

/* cosyrush.cc *****-C++-*****
**
**           Example program for the GAMMA Project
**
** The Program reads a spin system (number of spins, chemical shifts and
** J-coupling constants) from a file and simulates a normal COSY experiment
** using RuSH for quadrature detection in t1
**
** NOTE: Set to work only on HOMONUCLEAR systems
**
***** */

#include <gamma.h>
//
//           Define Constants
const int t1pts = 512;           // Number of t1 points
const int t2pts = 512;           // Number of t2 points
//
//           Begin Program
main (int argc, char* argv[])
{
    cout << "Homonuclear COSY Simulation with RuSH Mode\n";
//
//           Read in the Spin System
String filename;                 // Name of spin system file
query_parameter(argc, argv, 1,   // Get filename from command
    "nSpin system filename? ", filename); // line or ask for it
spin_system sys;                 // Declare spin system sys
sys.read(filename);              // Read system from filename
//
//           Set Offsets and Spectral Widths
double offset = sys.center();    // Find approx. spectrum center
sys.offsetShifts(offset);        // Offset shifts so centered
double NyqF = sys.Nyquist(0, 1.4); // Approximate Nyquist frequency
double t2dt = 1.0/(2.0*NyqF);    // t2 time increment
double t1dt = t2dt;              // t1 time increment
//
//           Set Up Hamiltonian
char J;
query_parameter(argc, argv, 2,   // Weak or strong coupling
    "nWeak or strong coupling (w/s)?", J);
gen_op H;
if (J == 'w')                 // Set Hamiltonian for
    H = Hcs(sys) + HJw(sys);   // strong or weak coupling
else
    H = Hcs(sys) + HJ(sys);
}

gen_op Upx = lxypuls_U(sys,0.0,90.0); // Propagator for x pulse
gen_op Upy = lxypuls_U(sys,90.0,90.0); // Propagator for y pulse
gen_op Ud1 = prop(H, t1dt);           // Propagator t1 delay increment
gen_op D = Fm(sys);                   // Detector to F-
gen_op sigma0, sigma1, sigma2, sigma3; // Set up density matrices
block_1D t2BLK(t2pts);                // Set 1D block for output
//
//           Pulse Sequence and I/O Setup
File cosyRe, cosyIm;                 // Declare and open two files
cosyRe.open("cosyRuSH.Redat",io_writeonly, a_create);
cosyIm.open("cosyRuSH.Imdat",io_writeonly, a_create);
//
//           Apply Pulse Sequence
gen_op sigma0 = sigma_eq(sys);        // Set density matrix equilibrium
gen_op sigma1 = evolve(sigma0, Upx);   // Apply first (PI/2)x pulse
gen_op sigma2 = sigma1;               // Initial sigma2 (t1 = 0)
for (int t1=0; t1<t1pts; t1++)        // Loop over all t1 increments
{
    sigma3 = evolve(sigma2, Upx);      // Apply second (PI/2)x pulse
    FID(sigma3,D,H,t2dt,t2pts,t2BLK); // Acquire data for reals
    Felix(cosyRe, t2BLK);              // Output block for reals: Felix
    sigma3 = evolve(sigma2, Upy);      // Apply second (PI/2)y pulse
    FID(sigma3,D,H,t2dt,t2pts,t2BLK); // Acquire data for imaginaries
    Felix(cosyIm, t2BLK);              // Output block for imgs: Felix
    evolve_ip(sigma2, Ud1);            // evolution to next t1
}
cosyRe.close();                      // Close files
cosyIm.close();
}

```

1.5.3 Discussion

Define System and NMR Parameters: The first two lines set the number of points to simulation along each dimension. These can be altered before compilation to suit individual needs. The program begins at the main statement, in this case main itself has arguments which allow the program to take use commands given on the command line. Following a brief output about the type of simulation, the spin system is read in from an external file. First a string variable “filename” is declared. The filename is set to be the first parameter given on the command line or, if none is given there, the user is asked for the filename. Having the spin system filename, a spin system variable “sys” is declared and read in from the file. Since the program is “automated” it finds an appropriate frequency offset (for better digital resolution) and shifts the spin system by this amount. An approximate Nyquist frequency is obtained from the function “Nyquist” and this is used to set the dwell times on both axes, “t1dt” and “t2dt” respectively.

Set up the Necessary Operators: It is left to the user whether the simulation should be done with strong or weak coupling. Again, this will be set from the second parameter on the command line or, if not given on the command line, the user will be prompted to input which is preferred. A general operator H is declared and set to be the sum of the isotropic shift Hamiltonian plus the chosen isotropic scalar coupling Hamiltonian. After this, three propagators (general operators) are created. These are Upx, Upy, and Ud1 for the 90x pulse, the 90y pulse, and the t_1 time increment - all of these are constant throughout the simulation and need to be computed only once. Also computed once is the general operator D used for the detection. The detector, D, is set to $F_- = F_x - iF_y$ so as to produce a complex data set along the t_2 axis. We then declare the four density matrices (corresponding to the pulse sequence diagram) and a data block “t2BLK” in which to collect the FID’s.

Apply the Pulse Sequence: The program sets the initial density matrix, σ_0 , to the equilibrium density matrix. A ninety degree pulse is then applied along the x-axis to produce σ_1 , and then the first σ_2 , $\sigma_2(t_1=0) = \sigma_1$, is set. Then begins the looping over all t_1 points desired. The density matrix $\sigma_3(t_1, x)$ is generated from applying a 90x pulse to $\sigma_2(t_1)$ and in turn $\sigma_3(t_1, x)$ is propagated in t_2 to produce $\sigma_3(t_1, t_2, x)$ and a corresponding FID which is output to the file cosyRe, “recosy.dat” in Felix compatible format. Then, the same $\sigma_2(t_1)$ is used to produce another σ_3 but this time with the application of a 90y pulse. Again $\sigma_3(t_1, y)$ is used to generate $\sigma_3(t_1, t_2, y)$ and a corresponding FID which is now output to the file cosyIm, “imcosy.dat”. The loop is then repeated for each t_1 but, rather than produce $\sigma_2(t_1)$ from σ_1 by evolving it for time t_1 we compute $\sigma_2(t_1)$ from $\sigma_2(t_1 - \Delta t_1)$ by evolving it for time Δt_1 . This is performed in the last step of the loop. The last two lines of the program close the two output files.

1.5.4 Felix Processing

The current version of Felix comes with a manual loaded with examples of macros for performing data processing. Unfortunately, most of these simply do not work at all¹. In this particular case, unlike processing a data set acquired using TPPI, we must take a complex Fourier transform along the t_1 axis. Since we also want to produce a contour plot of the processed real data we would prefer to maintain all data in a Felix real matrix. The Felix manual is very misleading and vague concerning this type of treatment so here we cover the required processing in fine detail. The following figure diagrams how the two data sets $S1$ and $S2$ are initially worked up in the t_2/F_2 dimension.

RuSH Processing in t_2/F_2 via Felix Real Matrix




Figure 0-9 Processing RUSH data sets in the T_2 dimension with Felix.

1. The manual currently being referred to is dated October 1990 and is for Felix Version 1.1. This manual exemplifies how extremely difficult it is to keep user documentation for a program current (hopefully GAMMA docs. aren't too similar in this point). We therefore tend to always include the Felix macros used in processing our simulated data in discussing these examples. From past experience, there is absolutely no guarantee that any of our listed macros will work in future Felix versions!

As each FID block is read in from the respective data file it is initially (**A**) zero filled to the Felix matrix dimension and then (**B**) Fourier transformed. The block is now complex frequency data and the reals and imaginaries contain redundant information since they differ only by a phase change. The block is also too large to fit into a row of the Felix matrix- the block dimension is correct but the matrix only has room for real data, not complex data. The block will be cut in half and fit correctly but must first be manipulated to isolate the useful data. In the next step (**C**), the reals are completely separated from the imaginaries. In the following step (**D**) the imaginaries are thrown out, they contain no new information, and the block not contains only real frequency information. The final step (**E**) puts the real frequency data into the Felix matrix row. Here, the blocks from the two data sets must be handled slightly differently (**E & E'**). We are storing all data in a real Felix matrix, yet we wish to perform a complex Fourier transform along t_1 . Since Felix always believes complex points are stored in (real, imaginary) pairs we must store our blocks in “real” block, “imaginary” block pairs. Since our “real” data comes from the *S1* data set and the “imaginary” data comes from the *S2* data set we must store the blocks in an interleaved order.

RuSH Processing in $t1/F1$ via Felix Real Matrix




Figure 0-10 Processing RUSH data sets in the $T1$ dimension with Felix.

We now discuss how Felix processes the data along the $t1/F1$ dimension. Each column is successively pulled out of the Felix real matrix (**A**). Since the matrix is real, the data will be real to Felix even though we have constructed the columns so they are in the Felix complex column format. The program must at this point be told the data is complex after which the data size will register as $1/2$ the value when removed from the matrix. The now complex data is zero filled to the matrix column

dimension (**B**) and Fourier transformed (**C**). As in the t2/F2 processing, we now have complex data which is too large to fit into the matrix but we only want the reals. So, we sort the reals and imaginaries (**D**) then throw away the imaginary part (**E**). Felix is told that the data is again real and it is then stored back into the matrix (**F**).

A conceptually simpler Felix workup can be performed using a complex array. The trade-off in simplicity is that the disk space used by the matrix twice that used in the macro manipulating the data in a real matrix. The diagram below shows how processing is done with the complex matrix.

RuSH Processing in t1/F1 via Felix Complex Matrix




Figure 0-11 Processing RUSH data sets in the T1 dimension with Felix.

Rows are successively read from each data set, $S1$ and $S2$ and Fourier transformed along the t_2 dimension (**A**). The real frequency data from $S1$ is blended with the imaginary frequency data from $S2$ (**B**) to form a new complex data set having the appropriate phase information for quadrature processing in t_1 . This is stored into the Felix complex matrix row (**C**). Afterwards, successive complex columns are read from the matrix (**D**), Fourier transformed and restored into the complex array (**E**).

1.5.5 Spin System

The spin system file, cosy.3.sys, is frequency shifted from the one used in the previous simulation. It is reproduced below.

SysName	(2) : COSYRuSH	- Name of the Spin System
NSpins	(0) : 3	- Number of Spins in the System
Iso(0)	(2) : 1H	- Spin Isotope Type
Iso(1)	(2) : 1H	- Spin Isotope Type
Iso(2)	(2) : 1H	- Spin Isotope Type
v(0)	(1) : 200.0	- Chemical Shifts in Hz
v(1)	(1) : 100.0	- Chemical Shifts in Hz
v(2)	(1) : 40.0	- Chemical Shifts in Hz
J(0,1)	(1) : 10.0	- Coupling Constants in Hz
J(0,2)	(1) : 0.0	- Coupling Constants in Hz
J(1,2)	(1) : 15.0	- Coupling Constants in Hz
Omega	(1) : 400	- Spectrometer Frequency in MHz (1H based)

1.5.6 Workup

We provide two Felix macros for data processing. As discussed in a previous section, we wish to use real matrices when working up the data and the corresponding macro is provided. A second macro is also given, one in which the data is stored in a complex Felix matrix. The former saves on disk space while the latter is perhaps more straight forward.

COSY RuSH Macro Using Felix Real Matrix

```
def redat cosyre                                ! set up a variable for the "real" data file
def imdat cosyim                                ! set up a variable for the "imag" data file
def matfile cosy                                ! set up a variable for the real Felix matrix file
def t1max 512
lb 10
ph0 0
ty Opening matrix file "&matfile " .....      ! report of macro progress
cmx                                              ! close any open matrix
cl                                              ! close any open files
mat &matfile write                              ! open the matrix file for writing
ty Using real data file "&redat " .....         ! another status report
ty Complex-FT of t2 into real F2 .....
for row 1 &t1max                                ! begin looping through the t1 values
re &redat                                       ! read in a complex FID from the "real" data
zf 1024
em                                              ! apodize the FID
ft                                              ! Fourier transform the complex data
ph
rev                                              ! reverse the frequencies for plotting
sep                                             ! split up all real/imaginary pairs
dat 0                                           ! call the data real, now twice a big
si 1024                                         ! use only the first half, all the true reals
sto 0 &nrow                                    ! store real row into real matrix - 1,3,5,7,9,...
ty reals row &row - matrix row &nrow
eva nrow (&nrow+2)                             ! increment matrix storage row count by two
next                                           ! back for next row of "real" input data
ty Using imaginary data file "&imdat " .
ty Complex-FT of t2 into imaginary F2 ...
def nrow 2
for row 1 &t1max                                ! reset the matrix storage row to be initially 2
re &imdat                                       ! read in a complex FID from the "imag" data
zf 1024
em                                              ! apodize the FID
ft                                              ! Fourier transform the complex data
rev                                             ! reverse the frequencies for plotting
exc
sep                                             ! split up all real/imaginary pairs
dat 0                                           ! call the data real, now twice a big
si 1024                                         ! use only the first half, all the true reals
sto 0 &nrow                                    ! store real row into real matrix - 2,4,6,8,...
ty imgs row &row - matrix row &nrow
eva nrow (&nrow+2)                             ! increment matrix storage row count by two
next                                           ! back for next row of "imag" input data
ty Complex-FT of t1 into F1 .....
```

```
for col 1 1024
  loa &col 0
  dat 1
  si 512
  zf 1024
  em
  ft
  ph
  rev
  sep
  dat 0
  si 1024
  sto &col 0
  ty col &col
  next
  ty Setting Up Contour Plot.....
  cmx
  mat &matfile read
  lim 1 1 1024
  lim 2 1 1024
  lvl 1e-4
  cpn 1
  pen 1
  nl 4
  cyc 2
  rmx 1 400 239 3 512 0 F2
  rmx 2 400 238 3 512 0 F1
  ty Plotting contours .....
  cp
  ty Processing finished .....
end
```

! load column of real data from real matrix
! but it is actually complex data, half the size

! apodize the FID
! Fourier transform the complex data

! reverse the frequencies for plotting
! split up all real/imaginary pairs
! again call the data real, now twice a big

! store the real data back into the real matrix

COSY RuSH Macro Using Felix Complex Matrix

```
def redat cosyre
def imdat cosyim
def matfile cosy
def t1max 512
def t2max 512
def total 512
lb 10
ph0 0
ty Opening matrix file "&matfile " .....
cmx
cl
mat &matfile write
ty Using real data file "&redat " .....
ty Complex-FT of t2 into real F2 .....
for row 1 &t1max
  re &redat
  em
  ft
  zi
  sto 0 &row
  ty reals row &row
  next
```

! set up a variable for the "real" data file
! set up a variable for the "imag" data file
! set variable for the complex Felix matrix file

! report of macro progress
! close any open matrix
! close any open files
! open the matrix file for writing
! another status report

! begin looping through the t1 values
! read in an FID from the "real" data
! apodize the FID
! Fourier transform the complex data
! zero the imaginaries, they'll be replaced
! store the data into the complex matrix
! tell which row is being processed
! go to the next row

```
ty Using imaginary data file "&imdat " .      ! another status report
ty Complex-FT of t2 into imaginary F2 ...      ! begin looping through the t1 values
for row 1 &t1max                                ! read corresponding matrix row or "reals"
loa 0 &row                                      ! store the complex row with reals in a buffer
stb 1                                          ! read corresponding row of "imags" data
re &imdat
si &t2max
em
ft
zr
adb 1                                          ! apodize the FID
ldb 1                                          ! Fourier transform the complex data
rev                                           ! zero the reals this time
sto 0 &row                                     ! add these imaginaries to the reals in buffer
ty imgs row &row                             ! get the new complex row from the buffer
next                                          ! reverse the frequency order for plotting
ty Complex-FT of t1 into F1 .....            ! store the row back into the matrix
for col 1 &total                             ! tell which row is being processed
loa &col 0                                    ! proceed to next row
si &t1max                                     ! report that the columns are ready
em                                           ! begin looping through the columns
ft                                           ! read in a complex column
ph 0                                          ! apodize the FID
rev                                           ! Fourier transform the complex data
sto &col 0                                   ! reverse the frequency order for plotting
ty col &col                                  ! store the complex column back in the matrix
next                                          ! report on the column progress
ty Setting Up Contour Plot.....              ! go get the next column
cmx                                           ! status report that plotting is ready
mat &matfile read                            ! first close the matrix
lim 1 1 512                                  ! reopen the matrix, now read only for safety!
lim 2 1 512                                  ! set the plotting limits
lvl 1e-4                                     ! set the initial contour intensity
cpn 1                                         ! contour both positive and negative peaks
pen 1                                         ! begin with pen 1 on first contour
nl 4                                          ! plot four contour levels
cyc 2                                         ! cycle the pen, + one color, - another color
rmx 1 300 1596 3 256 5.3 F2                 ! set the axes up
rmx 2 300 1596 3 256 5.3 F1
ty Plotting contours .....                  ! status report that plotting is starting
cp                                           ! draw a contour plot
ty Processing finished .....                ! status report that the macro is through
end                                           ! end the macro
```

1.5.7 Results

COSY RuSH Macro Using Felix Real Matrix

Once Felix has been started, the first task is to build a real matrix into which the data is to be placed. The command for this is (in Felix)¹

```
bld cosy 2 1024 1024 0
```

This creates a real matrix file called cosy.mat which is 1Kx1K. Now the macro listed previously can be run which will read the data in and process it. The command for this (still in Felix) is

```
ex recosy4
```

where recosy4.mac is the name of the file containing the macro². Once this macro has been executed in Felix, the COSY spectrum will automatically be displayed. To generate the hpgl contour plot file, the Felix commands are (once cp produces the screen plot)

```
hdv felix.hpgl      ! hardcopy device is file felix.hpgl
hpm 32              ! hardcopy plot mode is hpgl
hcp                 ! produce the plot
```

COSY RuSH Macro Using Felix Complex Matrix

Once Felix has been started, the first task is to build a complex matrix into which the data is to be placed. The command for this is (in Felix)³

```
bld cosy 2 512 512 1
```

This creates a complex matrix file called cosy.mat which is 512x512. Now the macro listed previously can be run which will read the data in and process it. The command for this (still in Felix) is

```
ex cmcosy4
```

where cmcosy4.mac is the name of the file containing the macro⁴. Once this macro has been executed in Felix, the COSY spectrum will automatically be displayed. To generate the hpgl contour plot file, the Felix commands are (once cp produces the screen plot)

```
hdv felix.hpgl      ! hardcopy device is file felix.hpgl
hpm 32              ! hardcopy plot mode is hpgl
hcp                 ! produce the plot
```

-
1. Currently the Felix documentation is contradictory in its use of bld. Experience says that a 0 at the end of the build statement produces a real array, a 1 or 2 produces a complex array.
 2. Felix may have trouble finding this macro file unless the "pre" command is used to specify which directory the macro files are located in. Keep in mind that Felix has problems (at least in UNIX) with capital letters in filenames and that directory names need to end with a /.
 3. See footnote 1.
 4. See footnote 4.

1.6 COSY with real pulses

1.6.1 Description

This example again performs a simple COSY simulation. In this instance, a TPPI¹ scheme is implemented in order to achieve quadrature detection along the t_1 axis. This is done by phase incrementing the initial pulse by 90 degrees on successive points along t_1 . For clarity, the pulse sequence diagram is repeated with density matrices used to coincide with variable names in the program code.




Figure 0-12 A simple COSY pulse sequence. Density operators label key points in simulation of this experiment.

In this example, rather than incrementing the receiver phase (normal TPPI), for convenience the first pulse phase is incremented. No relaxation effects will be considered and the pulses taken as ideal. The resulting output file is produced in Felix format and subsequently processed with that program.

1. The TPPI detection scheme is discussed in Derome, *Modern NMR Techniques for Chemistry Research*, page 83.

1.6.2 Program

```
/*cosy.2.cc*****-c++*-  
**  
**          COSY Example Using Rectangular Pulses          **  
**  
** Simulates a cosy experiment using rectangular (soft) pulses. No phase **  
** cycling is employed. The data is processed using NMRi from Tripos. **  
**  
*****/  
  
#include <gamma.h>  
main()  
{  
  const double dt1 = 0.001;          // increments in  $\tau_1$   
  const double dt2 = 0.001;          // increments in  $\tau_2$   
  const int pt1 = 1024;               // acquisition points for  $\tau_1$   
  const int pt2 = 1024;               // acquisition points for  $\tau_2$   
  const double puls_duration = 10.0e-6;  
  spin_system sys;                   // define the system;  
  gen_op sigma, sigma1, H, detect;  
  block_1D tmp(pt2);                 // variable for the response  
  block_2D data(pt1,pt2);            // variable for the 2D-data  
  sys.read("cosy.sys");               // read the system  
  sigma = sigma_eq(sys);              // equilibrium density matrix  
  H = Hcs(sys) + HJ(sys);             // Hamiltonian of the system)  
  detect = Fm(sys);  
  sigma = Sypuls(sys, sigma, H,"1H",  
                  0, puls_duration, 90); // apply the rf puls,0 Hz offset  
  for (int step=0; step<pt1; step++) // loop for  $\tau_1$   
  {  
    sigma1=evolve(sigma, H, step*dt1); // evolution during t1  
    sigma = Sypuls(sys, sigma, H,"1H",  
                  0, puls_duration, 90); // the second pulse  
    FID(sigma1, detect, H, dt2, pt2, tmp); // record the FID  
    data.put_block(step, 0, tmp);         // store it in the 2D block  
  }  
  NMRi ("cosy.dat", data);            // Output for the NMRi program  
}
```

1.6.3 Results

The data was processed with the NMRi program. The following processing was done:

- EI Set imaginary part to 0
- SB Phased Sine Bell (0.5, 1, 6)
- ZF Zero fill (No)
- FT Fourier transformation
- TP Transpose
- EI Set imaginary part to 0
- SB Phased Sine Bell (0.5, 1, 6)
 - ZFZero fill (No)
 - FTFourier transformation

The result is shown at the end of the chapter (Spectrum 3.1.2)

E-COSY Cross Peaks for Four Coupled Spins




Figure 0-14 The E-COSY cross peaks expected from four mutually coupled spins can be rationalized from the cross peaks in the three spin network. The eight peaks in the three spin system propagate into 2 sets of eight peaks which are shifted by the additional coupling constants to the two “active” spins A and B, namely J_{BD} and J_{AD} .

The pulse sequence and phase cycle for the E-COSY experiment are given in the following figure¹.

E-COSY Pulse Sequence and Phase Cycle




Figure 0-15 : Homonuclear E-COSY pulse sequence. The phase cycle is given in the previous JMR article in TABLE 1, page 477, with $N=K=3$.

This E-COSY sequence is identical to that used for MQF-COSY experiments. In MQF-COSY the phase cycle is adjusted to select out the coherence order desired. The E-COSY experiment is simply a linear combination of MQF-COSY so there is an additionally “weighting factor” W listed in the phase cycle.

1. A alternative E-COSY sequence is given in the previous *J. Chem. Phys* reference, page 6839, FIG. 4. That sequence has only two pulses, both of constant phase with the second pulse having a cycled pulse angle.

2.2 E-COSY Examples

The following set of example programs are covered in this chapter. They simulate E-COSY experiments and/or E-COSY related experiments and values. Both homonuclear and heteronuclear spin systems are treated.

Table 2: E-COSY Example Programs

Example	Page	Pulse ^a	Relaxation	System	Workup ^b
1. E-COSY	41	Ideal	No	Homonuclear	Felix
2. Alternate E-COSY	46	Ideal	No	Homonuclear	Felix
3. Complimentary	54	Ideal	No	Homonuclear	Felix
4. Heteronuclear E-COSY	61	Ideal	No	Heteronuclear	Felix
5. E-COSY Phase Cycles	69	N/A	N/A	N/A	FrameMaker
6. E-COSY with Relaxation	76	Ideal	Yes	Homonuclear	Felix
7. E-COSY with Relaxation		Ideal	Yes	Heteronuclear	Felix

a. Application of an “Ideal” pulse is “infinitely” short so no relaxation effects can be considered in such steps.

b. Workups labeled “Open” allow the user to choose the output type while the program is running. The output types currently available in GAMMA are FrameMaker, Felix, NMRi, and MATLAB

Each section in this chapter discusses at length a corresponding simulation in the above table. Although this document should contain listings of all files necessary to run the program, they should be included (reside on disk) with the GAMMA installation as well. These will lie in the */gamma/Example subdirectories.

2.3 E-COSY with Superoperators

2.3.1 Description

Our initial E- COSY simulation will implement the pulse sequence and phase cycle previously described. It is repeated below for reader clarity.

E-COSY Pulse Sequence and Phase Cycle





Figure 0-16 The standard E-COSY pulse sequence with appropriate phase cycle and weighting factors. The phase angles are given by $\phi = j(\pi/N)$ with $j = [0, 2N - 1]$; we have used $N=3$. The weighting factors are given by $W_0 = 4$ for $j=0$, $W_3 = 0$ for $j=K=3$, and $W_j = \frac{3}{4} \frac{(-1)^j}{\sin^2(\beta_j/2)}$ for all other j values. All formulae are obtained from the previous reference in JCP, page 6842. The weights are obtained with $B_0=B_1=0$ and the W_j are scaled by a factor of two.

In the figure, W is a weighting factor, ϕ the phase angle of the first two pulses, and ψ the phase angle of the detector. The evolution of the density matrix for this E-COSY sequence is depicted below in terms of Hilbert space propagators.

E-COSY Sequence in Terms of Hilbert Space Propagators



Accordingly, the evolution of the density matrix during an E-COSY experiment is given in the fol-

lowing set of equations.

$$\sigma_1(\phi) = U_p(x, \phi) \sigma_{eq} U_p^{-1}(x, \phi) = R_z(\phi) U_p(x) R_z^{-1}(\phi) \sigma_{eq} R_z(\phi) U_p^{-1}(x) R_z^{-1}(\phi) \quad (\text{EQ 1})$$

$$\sigma_2(t_1, \phi) = U_d(t_1) \sigma_1(\phi) U_d^{-1}(t_1) \quad (\text{EQ 2})$$

$$\begin{aligned} \sigma_3(t_1, \phi) &= U_p(-x) U_p(x, \phi) \sigma_2(t_1, \phi) U_p^{-1}(x, \phi) U_p^{-1}(-x) \\ &= U_p(-x) R_z(\phi) U_p(x) R_z^{-1}(\phi) \sigma_2(t_1, \phi) R_z(\phi) U_p^{-1}(x) R_z^{-1}(\phi) U_p^{-1}(-x) \end{aligned} \quad (\text{EQ 3})$$

$$\sigma(t_1, t_2, \phi) = U_d(t_2) \sigma_3(t_1, \phi) U_d^{-1}(t_2) \quad (\text{EQ 4})$$

The FID is then computed by performing trace operations with a detection operator, F_- , which is also phase cycled.

$$FID(t_1, t_2, \phi, \psi) = \text{Tr}\{F_-(\psi) \sigma(t_1, t_2, \phi)\} \quad (4-1)$$

The last step is that we sum the FID's over the phase cycle

$$FID(t_1, t_2) = \sum_{\phi, \psi}^{\text{cycle}} FID(t_1, t_2, \phi, \psi) \quad (4-2)$$

With GAMMA, we could proceed implementing these equations as outlined and directly perform an E-COSY simulation. Such a task would require some scheme for performing the phase cycle, perhaps involving a complicated algorithm and definitely involving repetitive computations. Thus, rather than using these six equations as written, we will shuffle around the formulae and attempt to isolate components involving the phase cycle from components involving the t_1 and t_2 time increments. The result will simplify implementation of the E-COSY sequence. From the first two equations we have

$$\sigma_2(t_1, \phi) = U_d(t_1) R_z(\phi) U_p(x) R_z^{-1}(\phi) \sigma_{eq} R_z(\phi) U_p^{-1}(x) R_z^{-1}(\phi) U_d^{-1}(t_1) \quad .$$

Rotations about the z-axis cannot affect the equilibrium density matrix, so the inner rotation operators may be safely removed. Furthermore, rotations about z commute with the time evolution propagators due to $[H_0, F_z] = 0$, so we may switch the order of these operations. The result is

$$\sigma_2(t_1, \phi) = R_z(\phi) U_d(t_1) U_p(x) \sigma_{eq} U_p^{-1}(x) U_d^{-1}(t_1) R_z^{-1}(\phi) \quad . \quad (4-3)$$

We define a new density matrix operators which are independent of the phase angle ϕ ,

$$\sigma_1 = U_p(x) \sigma_{eq} U_p^{-1}(x) \quad (4-4)$$

and

$$\sigma_2(t_1) = U_d(t_1) \sigma_1 U_d^{-1}(t_1) \quad , \quad (4-5)$$

so that the previous formula for $\sigma_2(t_1, \phi)$, (4-3), becomes simply

$$\sigma_2(t_1, \phi) = R_z(\phi)\sigma_2(t_1)R_z^{-1}(\phi) \quad . \quad (4-6)$$

Keep in mind our objectives: to separate time evolution steps from phase cycle steps. Clearly, equation (4-6) shows how to obtain $\sigma_2(t_1, \phi)$ from independent phase and time evolutions steps.

A cancellation of the rotation operators in (4-6) occurs when the next step is formulated. From (EQ 3),

$$\begin{aligned} \sigma_3(t_1, \phi) &= U_p(-x)R_z(\phi)U_p(x)R_z^{-1}(\phi)\sigma_2(t_1, \phi)R_z(\phi)U_p^{-1}(x)R_z^{-1}(\phi)U_p^{-1}(-x) \\ \sigma_3(t_1, \phi) &= U_p(-x)R_z(\phi)U_p(x)\sigma_2(t_1)U_p^{-1}(x)R_z^{-1}(\phi)U_p^{-1}(-x) \end{aligned} \quad (4-7)$$

Again, note that the phase cycle is explicitly removed from the t_1 incrementation at this point and we are at the last step in the pulse sequence, detection.

Detection involves propagation for the t_2 delay followed by multiplication with a detection operator respectively. We have

$$F_-(\psi)\sigma(t_1, t_2, \phi) = R_z(\psi)F_-R_z^{-1}(\psi)\sigma(t_1, t_2, \phi) \quad . \quad (4-8)$$

and can utilize the relationship $Tr\{AB\} = Tr\{BA\}$, to produce

$$\begin{aligned} Tr\{F_-(\psi)\sigma(t_1, t_2, \phi)\} &= Tr\{R_z(\psi)F_-R_z^{-1}(\psi)\sigma(t_1, t_2, \phi)\} \\ &= Tr\{F_-R_z^{-1}(\psi)\sigma(t_1, t_2, \phi)R_z(\psi)\} \end{aligned}$$

If the delay t_2 is explicitly written we obtain

$$\begin{aligned} Tr\{F_-(\psi)\sigma(t_1, t_2, \phi)\} &= Tr\{F_-R_z^{-1}(\psi)U_d(t_2)\sigma_3(t_1, \phi)U_d^{-1}(t_2)R_z(\psi)\} \\ &= Tr\{F_-U_d(t_2)R_z^{-1}(\psi)\sigma_3(t_1, \phi)R_z(\psi)U_d^{-1}(t_2)\} \end{aligned} \quad (4-9)$$

Thus, the phase cycle over ψ has been isolated from the detection and t_2 incrementation. If we again focus on the components which depend upon the phase cycle, we have

$$\begin{aligned} R_z^{-1}(\psi)\sigma_3(t_1, \phi)R_z(\psi) &= \sigma_3(t_1, \phi, \psi) \\ &= R_z^{-1}(\psi)U_p(-x)R_z(\phi)U_p(x)\sigma_2(t_1)U_p^{-1}(x)R_z^{-1}(\phi)U_p^{-1}(-x)R_z(\psi) \end{aligned} \quad . \quad (4-10)$$

The entire E-COSY pulse sequence phase cycle is contained in the above equation and we have fulfilled our objective of completely separating the time incrementation from the phase cycle mathematically. We can formulate a new propagator which contains the phase cycle. Letting

$$U_{mix}(\phi, \psi) = R_z^{-1}(\psi)U_p(-x)R_z(\phi)U_p(x) \quad , \quad (4-11)$$

the trace equation becomes,

$$\begin{aligned} Tr\{F_-(\psi)\sigma(t_1, t_2, \phi)\} &= Tr\{F_-U_d(t_2)U_{mix}(\phi, \psi)\sigma_2(t_1)U_{mix}^{-1}(\phi, \psi)U_d^{-1}(t_2)\} \\ &= Tr\{F_-U_d(t_2)\sigma_3(t_1, \phi, \psi)U_d^{-1}(t_2)\} \end{aligned} \quad (4-12)$$

To clarify what has been accomplished, the E-COSY pulse sequence (for mathematical implementation, not experimental implementation) is rewritten in the following diagram to show the distinction between the time delays and the phase cycle steps.

Equivalent E-COSY Sequence in Terms of Hilbert Space Propagators




Figure 0-17 The E-COSY pulse sequence in terms of adjusted Hilbert space propagators. At this point the propagators for the delays are separate from the propagators involving the phase angles.

Of course, in order to compute our 2-dimensional E-COSY spectrum we need to sum FID's over each phase cycle according to (4-2)

$$FID(t_1, t_2) = \sum_{\phi, \psi}^{\text{cycle}} FID(t_1, t_2, \phi, \psi)$$

and we have a different mixing propagator for each combination of ϕ and ψ .

We now employ the power of superoperators and it will hopefully become clear to the reader why we have bothered rewriting the mathematical steps describing E-COSY. We shall form the unitary transformation superoperator equivalent to the mixing propagator. The equation then becomes,

$$Tr\{F_-(\psi)\sigma(t_1, t_2, \phi)\} = Tr\left\{F_-U_d(t_2)[\hat{\Gamma}_{mix}(\phi, \psi)\sigma_2(t_1)]U_d^{-1}(t_2)\right\} \quad (4-13)$$

and we can take advantage of superoperator linearity in steps which are not linear in a Hilbert space

propagator formulation. Now when we sum over the phase cycle we have

$$\begin{aligned}
 FID(t_1, t_2) &= \sum_{\phi, \psi} FID(t_1, t_2, \phi, \psi) \\
 FID(t_1, t_2) &= \sum_{\phi, \psi}^{cycles} Tr \left\{ F_- U_d(t_2) [\hat{\Gamma}(\phi, \psi) \sigma_2(t_1)] U_d^{-1}(t_2) \right\} \\
 FID(t_1, t_2) &= Tr \left\{ F_- U_d(t_2) \left[\sum_{\phi, \psi}^{cycles} \hat{\Gamma}(\phi, \psi) \sigma_2(t_1) \right] U_d^{-1}(t_2) \right\} \\
 FID(t_1, t_2) &= \left\{ F_- U_d(t_2) [\hat{\Gamma}_{cycle}(\phi, \psi) \sigma_2(t_1)] U_d^{-1}(t_2) \right\}
 \end{aligned} \tag{4-14}$$

The phase cycle needed for the E-COSY simulation is now entirely contained in the superoperator, and may be applied in a single step. On a computational level, this means that the phase cycling loop can be removed from the loop over t_1 increments. The code becomes more concise and the computation more efficient¹. The E-COSY pulse sequence in terms of the phase cycle superoperator can be depicted as follows.

E-COSY Sequence Using A Superoperator Phase Cycle





Figure 0-18 The E-COSY pulse sequence in terms of Hilbert space propagators and a superoperator for the phase cycle. The entire phase cycle which is essential for E-COSY has been replaced by a single evolution under a Liouville space superoperator.

The simulation will be performed exactly this in this manner. We initially form the equilibrium density matrix, the pulse propagators, the t_1 propagator for t_1 incrementation time, and the superoperator representing the entire E-COSY phase cycle.

1. This is true only for small spin systems, <5 spins. As the system size increases it becomes more difficult for computers to handle superoperators due to the size of the Liouville space.

Program Outline for the E-COSY Sequence

Step	Density Matrix
1. Form Equilibrium Density Matrix	σ_{eq}
2. Apply the 1st x-pulse	σ_1
3. Evolve for t_1 time	$\sigma_2(t_1)$
4. Evolve under E-COSY “Phase Cycle”	$\sigma_3(t_1)$
5. Acquire FID for this t_1 time	$\sigma(t_1, t_2)$
6. Repeat Steps 3-6 for each t_1 time	



2.3.2 Program

```

/* ecosy1.cc *****-C++-
**
**          Example Program for the GAMMA Library
**
** This program simulates an E-COSY experiment with complete phase
** cycling. The pulses are ideal and relaxation effects are not considered.
**
** See: C. Griesinger, O.W. Sorensen, and R.R. Ernst, JMR, 75, 474-492,
** (1987), "Practical Aspects of the E-COSY Technique. Measurement of
** Scalar Spin-Spin Coupling Constants in Peptides". The required phase
** cycle is found on page 477 of this reference.
**
** Note: This program treats only homonuclear spin systems.
** Superoperators are used to implement the phase cycle
** which speeds up the computation for small systems
** but this program is cumbersome for large systems.
**
***** */

#include "gamma.h"

//          Define Constants

const P_cycl = 12;          // Phase cycle length
const double P_mix[12] =   // Pulse phase cycle
    {0,60,0,60,0,60,0,300,120,300,240,300};

//          Begin Program

main (int argc, char* argv[])
{
    cout << "\n\tSimulation of an E-COSY spectrum in TPPI mode\n";
    cout << "\t\t (Hard Pulses, No Relaxation)\n";

    Read in Parameters

    String filename;          // Name of spin system file
    query_parameter(argc, argv, 1, // Get filename from command
        "\n\tSpin system filename? ", filename); // line or ask for it
    spin_system sys;          // Declare spin system sys
    sys.read(filename);        // Read system from filename
    int t1pts, t2pts;
    query_parameter(argc, argv, 2, // Get number FID of points
        "\n\tAcquisition Size? ", t2pts);
    t1pts = 2*t2pts;          // Set t1 size for TPPI
    String J;
    query_parameter(argc, argv, 3, // Weak or strong coupling
        "\n\tWeak or strong coupling (w/s)? ", J);
    double offset = query_offset(sys, 0, 1); // Ask for an offset frequency
    double NyqF = query_Nyquist(sys, 0);    // Choose a Nyquist frequency
    double t2dt = 1.0/(2.0*NyqF);           // Dwell time, quadrature
    double t1dt = t2dt/2.0;                // t1 time increment, TPPI

```

Set Up Operators, Superoperator

```

gen_op H;          // Set Hamiltonian for
if (J == "w") H = Hcs(sys) + HJw(sys) // Strong or weak coupling
else H = Hcs(sys) + HJ(sys);
gen_op Ud1 = Rz(sys,+90)*prop(H,t1dt); // Delay t1 + TPPI propagator
gen_op D = Fm(sys); // Detector to F-
gen_op Upx = lxypuls_U(sys,0.,90.); // Propagator for x pulse
gen_op Upmx = lxypuls_U(sys,0.,-90.); // Propagator for -x pulse
gen_op U_mix; // Temporary mixing propagator
super_op G_mix; // Phase cycle superoperator
gen_op sigma3;

//          Construct Mixing, Phase Cycle Superoperator

double conv = acos(-1)/180;
double P_det = 0;
cout << "\n\tE-COSY Phase Cycle\n";
cout << "\n\t# scan          betareference\n";
for ( int i=0; i<P_cycl; i++ )
{
    cout << "\t" << i << "\t" << P_mix[i] // Output E-COSY phase
        << "\t" << P_det << "\n";
    U_mix = Rz(sys,-P_det); // Detector phase cycle
    U_mix *= Upmx; // 3rd 90 Pulse (-x)
    U_mix = Rz(sys,+P_mix[i]); // Phase shift pulses 1 & 2
    U_mix *= Upx; // 2nd 90 Pulse
    U_mix.Op_base(H); // Put in eigenbasis of Ho
    G_mix += U_transform(U_mix); // Add to U transform superop
    P_det = acos(-cos(P_det*conv))/conv; // Adjust detector phase
}

//          Apply Pulse Sequence

File ecosy; // Declare and open file
ecosy.open("ecosy.dat",
    io_writeonly, a_create);
block_1D t2BLK(t2pts); // Set 1D block for output
gen_op sigma0 = sigma_eq(sys); // Set density matrix equilibrium
gen_op sigma1 = evolve(sigma0,Upx); // Apply first (PI/2)x pulse
gen_op sigma2 = sigma1; // Initial sigma2 (t1 = 0)
for (int t1=0; t1<t1pts; t1++) // Loop over all t1 increments
{
    sigma3 = evolve(sigma2, G_mix); //superop phase cycle
    FID(sigma3,D,H,t2dt,t2pts,t2BLK); //acquisition
    Felix(ecosy, t2BLK); //output block: Felix
    evolve_ip(sigma2, Ud1); //evolution next t1
}
ecosy.close(); // Close file
double Om = sys.Omega(); // Spectrometer frequency
Felix2D_params(cout, Om, 2.0*NyqF,
    2*t2pts, offset); // Output Felix parameters
cout << "\n\n"; // Keep screen nice
}

```

2.3.3 Discussion

Define Constants: In this program the E-COSY phase cycle are explicitly defined; A 12 step cycle and associated phases in accordance with Figure 0-16.

Read in Parameters: As in prior simulations the spin system is read in from a disk file, the name may be supplied directly when running the program. The user is also asked for the acquisition size, whether to use a strong or weak coupling Hamiltonian, any desired offset, and a Nyquist frequency. Note that because TPPI will be used for phase sensitivity along the t_1/f_1 axis, twice as many points are taken along this axis that along the t_2 axis, and the t_1 incrementation time is half of the dwell time.

Set up Operators, Superoperator: The Hamiltonian is first formulated for either strong or weak coupling depending upon the input choice. Next, the t_1 time increment propagator is formulated. Here, in order to produce an acquisition with TPPI, a 90 degree rotation about the z-axis is added so that each incremented step includes the TPPI phase shift. The detector is set to F- in the next line. Propagators for the two 90 pulses are now computed, the first along x and the second along -x. A temporary operator is declared for used in summing over the phase cycle. The next line declares a superoperator which will ultimately account for the phase cycle. Finally, a working density matrix is specified.

Construct Mixing, Phase Cycle Superoperator: As outlined mathematically at the start of this section, the entire phase cycle can be removed from the t_1 and t_2 looping by use of a superoperator. In this section, the phase cycle loop is applied and the appropriate superoperator constructed. The loop goes over the twelve steps of the cycle (defined in the constants section) and at each step the phases are output to the terminal. Initially the detector phase is set to zero. The propagators for each step are then multiplied together in the opposite order of the pulse sequence due to the ordering of the unary *= step. At the end of each cycle, the unitary transform superoperator is summed to account for the mixing propagator. The last step in this section adjusts the detector phase to +/- 180.

Apply Pulse Sequence: Since the superoperator intrinsically contains the phase cycling, the pulse sequence application is quite simple. A file name ecosy is created and opened. A 1D data block called t2BLK is constructed to temporarily store the FID. Following the pulse sequence diagram, the initial density matrix is set to equilibrium. The next step is to apply a 90 degree pulse (without phase cycling) along the x-axis. This is evolved next by the propagator which accounts for t_1 incrementation as well as TPPI phase adjustments, but it is not done on the first step through the t_1 increments, at $t_1=0$. Rather, this step is done at the end of the loop. The next step is evolution under the superoperator which contains all the phase cycling. This density matrix, sigma3, is then used to simulate an FID which is subsequently sent to the file ecosy.dat in Felix format. Again, the last step of the loop is actually the earlier t_1 incrementation. The last program steps close the file, get the spectrometer frequency, output Felix parameters necessary for spectral workup, and clear the screen.

2.3.4 Spin System

For this example we will first treat a simple 3 spin proton system in order to reproduce the published simulations¹ and verify that our program is producing the proper results.

Initial 3-Spin Proton System for E-COSY Verification




Figure 0-19 Simple 3 spin protons systems to be utilized in verification of the E-COSY simulation program.

ecosy1A.sys

SysName	(2) : ECOSY1A	- Name of the Spin System
NSpins	(0) : 3	- Number of Spins in the System
Iso(0)	(2) : 1H	- Spin Isotope Type
Iso(1)	(2) : 1H	- Spin Isotope Type
Iso(2)	(2) : 1H	- Spin Isotope Type
PPM(0)	(1) : -.1	- Chemical Shifts in PPM
PPM(1)	(1) : 0.0	- Chemical Shifts in PPM
PPM(2)	(1) : .1	- Chemical Shifts in PPM
J(0,1)	(1) : 8.0	- Coupling Constants in Hz
J(0,2)	(1) : 3.0	- Coupling Constants in Hz
J(1,2)	(1) : 2.5	- Coupling Constants in Hz
Omega	(1) : 500	- 1H Spectrometer Frequency in MHz

ecosy1B.sys = ecosy1A.sys with replacements below

J(1,2)	(1) : -.75	- Coupling Constants in Hz
--------	------------	----------------------------

Afterwards we shall examine a system of both 3 and 4 spins which are meant to represent a prototypical amino acid.

2.3.5 Workup

The 2D-data sets produced from this simulation were processed with the program Felix. Each data set was subjected to the following Felix macro, called ecosy1.mac (the !comments to the right are not to be included in the macro - these are no longer allowed in the current Felix version.). Note that this macro *assumes that the simulation had an acquisition size of 512 points*.

1. C. Griesinger, O.W. Sørensen, and R.R. Ernst (1987), *JMR*, **75**, 474-492, "Practical Aspects of the E.COSY Technique. Measurement of Scalar Spin-Spin Coupling Constants in Peptides". These systems are those used in this paper to generate their E-COSY spectra in FIG. 15., page 488.

```
def datfile ecosy      ! simulated data file name "ecosy.dat"
def matfile ecosy      ! matrix file name "ecosy.mat"
def t1max 1024         ! dimension t1 size
def t2max 512          ! dimension t2 size
def total 1024
ty Opening source file " &datfile " .....! comment showing
which data file
ty Opening matrix file " &matfile " .....! comment showing
which matrix file
cmx                   ! close any existing matrix
cl                   ! close any existing dat files
mat &matfile write    ! open matrix ecosy.mat
ty Filling matrix with data ...! comment that data processing
starting
lb 4                  ! set line broadening for apodization
ph0 90               ! set zero order phase correction
ph1 0                ! set first order phase correction
ty Complex-FT of t2 into F2 ...! comment starting F2 transfor-
mations
for row 1 &t1max      ! loop through all t1 rows
re &datfile          ! read block of ecosy.dat file
si &t2max             ! set the size
em                   ! exponential multiplication
sb &t2max 90          ! sine bell windowing function
zf &total            ! zero fill to total size
ft                   ! complex FFT
red                  ! reduce data to real
sto 0 &row            ! store processed data into matrix
next                 ! loop back for next block
ty Real-FT of t1 into F1 ! comment for t1 workup
for col 1 &total      ! loop through all columns
loa &col 0            ! retrieve column from matrix
si &t1max             ! specify the size
em                   ! exponential multiplication
sb &t1max 90          ! sine bell
zf 2048              ! zero fill (will loose with real fft)
rft                  ! real FFT
ph                   ! apply phase adjustment
red                  ! reduce data to real
sto &col 0            ! store back into matrix
next                 ! get next column
ty Plotting contours ..... ! comments for processing point
```

```
ty zooming cross peaks ..
cmx                   ! clear matrices
mat &matfile read     ! open ecosy.mat for reading only
lim 1 1 1024          ! set plot limits
lim 2 1 1024
lvl 5e-4              ! set scaling
cpn 1                 ! contour positive and negative peaks
pen 1                 ! set first pen to pen 1
nl 4                  ! four contour levels
cyc 2                 ! pen cycle
cp                    ! contour plot
ty Processing finished .. ! comment for processing end
end
```

Before running this macro, construct matrix ecosy.mat with cmd


bld ecosy 2 1024 1024 0! real array
which constructs a 2 dimensional array (1K x 1K) of real numbers.
After processing the matrix axes are set with the rmx command
(used twice, once for each axis).

```
rmx
1(2)                  ! specify the dimension (1 or 2)
500                   ! spectrometer frequency
128                   ! spectral width
3                     ! PPM plotted on axis
512                   ! reference point
0                     ! reference point value
F2(F1)                ! axis label (F2 or F1)
```

The simulated E-COSY spectra were placed into this FrameMaker document in encapsulated Postscript format which resulted from sending the HPGL output of Felix through the filter program hpgl-toeps provided by FrameMaker. To generate the hpgl contour plot file, the Felix commands are (once cp produces the screen plot)

```
hdv felix.hpgl        ! hardcopy device is
file felix.hpgl
hpm 32                 ! hardcopy plot mode
is hpgl
hpc                     ! produce the plot
```

2.3.6 Results

Simulated E-COSY Spectra on a 3-Spin Proton System

The simulation was repeated for the two 3-spin systems, the input spin systems files are used to label the plots. Each plot was output from Felix, placed into FrameMaker MIF format, rescaled and annotated, and finally placed onto this page. All were performed with weak coupling.

ecosy13D.sys = ecosy13A.sys with replacements below

J(0,2)	(1) : 4.9	- Coupling Constants in Hz
ecosy14A.sys		
SysName	(2) : system	- Name of the Spin System
NSpins	(0) : 4	- Number of Spins in the System
Iso(0)	(2) : 1H	- Spin Isotope Type
Iso(1)	(2) : 1H	- Spin Isotope Type
Iso(2)	(2) : 1H	- Spin Isotope Type
Iso(3)	(2) : 1H	- Spin Isotope Type
PPM(0)	(1) : -.420	- Chemical Shifts in PPM
PPM(1)	(1) : -.035	- Chemical Shifts in PPM
PPM(2)	(1) : .250	- Chemical Shifts in PPM
PPM(3)	(1) : .385	- Chemical Shifts in PPM
J(0,1)	(1) : 6.8	- Coupling Constants in Hz
J(0,2)	(1) : 0.0	- Coupling Constants in Hz
J(0,3)	(1) : 0.0	- Coupling Constants in Hz
J(1,2)	(1) : 9.8	- Coupling Constants in Hz
J(1,3)	(1) : 0.0	- Coupling Constants in Hz
J(2,3)	(1) : -15	- Coupling Constants in Hz
Omega	(1) : 256	- Spectrometer Frequency in MHz (1H based)

ecosy14B.sys = ecosy14A.sys with replacements below

J(0,3)	(1) : -.9	- Coupling Constants in Hz
---------------	------------------	-----------------------------------

ecosy14C.sys = ecosy14A.sys with replacements below

J(1,3)	(1) : 3.9	- Coupling Constants in Hz
---------------	------------------	-----------------------------------

ecosy14D.sys = ecosy14A.sys with replacements below

J(0,3)	(1) : -.9	- Coupling Constants in Hz
J(1,3)	(1) : 3.9	- Coupling Constants in Hz

2.3.8 Workup


The 2D-data sets were processed in the same manner as for the previous plots except that the data matrices were kept smaller. Due to the inherent limitations in the Felix macro used, a second macro is supplied, *ecosy1b.mac* which *assumes that the simulation had an acquisition size of 128 points*.

Prior to running this Felix macro the matrix *ecosy.mat* must be constructed with the following Felix command

```
bld ecosy 2 256 256 0          ! real array
```

2.3.9 Results

Simulated E-COSY Spectra on a 3,4-Spin “Amino Acid” Proton System



2.4 Alternative E-COSY

2.4.1 Description

An alternative pulse sequence has been given in the literature as being equivalent to the E-COSY sequence used in the last simulation¹. The pulse sequence is more simplistic and, although perhaps more difficult to implement experimentally, readily implemented in a GAMMA program.

Alternative E-COSY Pulse Sequence





Figure 0-21 The alternate E-COSY pulse sequence with appropriate phase cycle and weighting factors. For $\psi = 180$, the weights may be replaced by negative values (i.e. subtract FID's).

The evolution of the density matrix for this E-COSY sequence is depicted below in terms of Hilbert space propagators.

Alternative E-COSY Sequence in Terms of Hilbert Space Propagators



Density matrix evolution through this E-COSY sequence is given in the following set of equations.

1. C. Griesinger, O.W. Sørensen, and R.R. Ernst (1986), *J. Chem. Phys.*, **85**, 6837-6851, "Correlation of connected transitions by two-dimensional NMR spectroscopy". See FIG. 4 on page 6839.

$$\sigma_1 = U_p(x)\sigma_{eq}U_p^{-1}(x) \quad (\text{EQ 5})$$

$$\sigma_2(t_1) = U_d(t_1)\sigma_1U_d^{-1}(t_1) \quad (\text{EQ 6})$$

$$\sigma_3(t_1, \phi) = U_p(y, \phi)\sigma_2(t_1)U_p^{-1}(y, \phi) \quad (\text{EQ 7})$$

$$\sigma(t_1, t_2, \phi) = U_d(t_2)\sigma_3(t_1, \phi)U_d^{-1}(t_2) \quad (\text{EQ 8})$$

The FID is then computed by performing trace operations with a detection operator, F_- , which is also phase cycled.

$$FID(t_1, t_2, \phi, \psi) = Tr\{F_-(\psi)\sigma(t_1, t_2, \phi)\} \quad (8-1)$$

The last step is that we sum the FID's over the phase cycle

$$FID(t_1, t_2) = \sum_{\substack{\text{cycle} \\ \phi, \psi}} FID(t_1, t_2, \phi, \psi) \quad (8-2)$$

Unlike the previous example, the step implementing the phase cycle is virtually independent from the steps involving the t_1 and t_2 time increments. Use of superoperator will not be done but we will still isolate the phase cycle. The first mathematical adjustment will be the realization that the detection phase cycle simply adds (0) or subtracts (180) the FID, and that we can remove the psi angle from the detection operator.

$$FID(t_1, t_2, \phi, \psi) = Tr\{F_-(\psi)\sigma(t_1, t_2, \phi)\} = (-1)^{\psi/\pi} Tr\{F_-\sigma(t_1, t_2, \phi)\}.$$

Expanding this result based on the previous equations produces

$$FID(t_1, t_2, \phi, \psi) = (-1)^{\psi/\pi} Tr\left\{F_-U_d(t_2)\sigma_3(t_1, \phi)U_d^{-1}(t_2)\right\}.$$

$$FID(t_1, t_2, \phi, \psi) = (-1)^{\psi/\pi} Tr\{F_-U_d(t_2)U_p(y, \phi)\sigma_2(t_1)U_p^{-1}(y, \phi)U_d^{-1}(t_2)\}$$

For the E-COSY we must sum over the phase cycle, weighting each by an appropriate amount.

$$FID(t_1, t_2) = \sum_{\substack{\text{cycle} \\ \phi, \psi}} (-1)^{\psi/\pi} Tr\{F_-U_d(t_2)U_p(y, \phi)\sigma_2(t_1)U_p^{-1}(y, \phi)U_d^{-1}(t_2)\}$$

$$FID(t_1, t_2) = Tr\left\{F_-U_d(t_2) \sum_{\substack{\text{cycle} \\ \phi, \psi}} W_\phi (-1)^{\psi/\pi} [U_p(y, \phi)\sigma_2(t_1)U_p^{-1}(y, \phi)] U_d^{-1}(t_2)\right\} \quad (8-3)$$


The factor (-1) in this scheme can be included in the weighting factor and our equation simplifies to

$$FID(t_1, t_2) = Tr \left\{ F U_d(t_2) \sum_{\phi}^{cycle} W_{\phi} [U_p(y, \phi) \sigma_2(t_1) U_p^{-1}(y, \phi)] U_d^{-1}(t_2) \right\}, \quad (8-4)$$

and the necessary weights, W_{ϕ} , are the +/- those previously given; their sign depends upon the angle phi (or correspondingly the detector angle). The simulation will be performed in this manner.

Program Outline for the Alternative E-COSY Sequence

Step	Density Matrix
-------------	-----------------------

- | | |
|---|--------------------|
| 1. Form Equilibrium Density Matrix | σ_{eq} |
| 2. Apply the 1st 90 x-pulse | σ_1 |
| 3. Evolve for t_1 time | $\sigma_2(t_1)$ |
| 4. Apply phi pulse, E-COSY Phase Cycle | $\sigma_3(t_1)$ |
| 5. Acquire FID for this t_1 time | $\sigma(t_1, t_2)$ |
| 6. Repeat Steps 3-6 for each t_1 time | |
- 

Note that this scheme is equivalent to the outline used in the previous example. The difference lies in the implementation of Step 4. Previously, a superoperator was used in implementing the E-COSY phase cycle. Now, after t_1 evolution, the system will be evolved through y-pulses of differing angles and the results scaled by a weighting factor then summed.

2.4.2 Program

```
/* ecosy2.cc *****-C++-*****
**
**           Example Program for the GAMMA Library
**
** This program simulates either an E-COSY experiment with
** complete phase cycling. The pulses are taken to be ideal
** and the effects of relaxation are not considered.
**
** See: C. Griesinger, O.W. Sorensen, and R.R. Ernst, JCP,
** 85(12), 6837-6852, (1986), "Correlation of Connected
** Transitions by Two-Dimensional NMR Spectroscopy". The
** required pulse sequence is found on page 6839 of this
** reference.
**
** Note: This program treats only homonuclear spin systems.
**
***** */

#include "gamma.h"
//
//           Define Constants
//
const P_cycl = 5; // Phase cycle length
const double P_mix[5] = {0,60,120,240,300}; // Pulse phase cycle, E-COSY
const double P_det[5] = {4,-3,1,1,-3}; // Detector phase cycle, E-COSY
//
//           Begin Program
//
main (int argc, char* argv[])
{
    cout << "\n\tAlternate Simulation of E-COSY\n";
    cout << "\t (TPPI Mode, Hard Pulses, No Relaxation)\n";
    //           Read in Parameters
    String filename; // Name of spin system file
    query_parameter(argc, argv, 1, // Get filename from command
        "\n\tSpin system filename? ", filename); // line or ask for it
    spin_system sys; // Declare spin system sys
    sys.read(filename); // Read system from filename
    int t1pts, t2pts;
    query_parameter(argc, argv, 2, // Get number FID of points
        "\n\tAcquisition Size? ", t2pts);
    t1pts = 2*t2pts; // Set t1 size for TPPI
    String J;
    query_parameter(argc, argv, 3, // Weak or strong coupling
        "\n\tWeak or strong coupling (w/s)? ", J);
    double offset = query_offset(sys, 0, 1); // Ask for an offset frequency
    double NyqF = query_Nyquist(sys, 0); // Choose a Nyquist frequency
    double t2dt = 1.0/(2.0*NyqF); // Dwell time, quadrature
    double t1dt = t2dt/2.0; // t1 time increment, TPPI
}
```

Set Up Operators, Superoperator

```
gen_op H; // Set Hamiltonian for
if(J == "w") // Strong or weak coupling
    H = Hcs(sys) + HJw(sys);
else
    H = Hcs(sys) + HJ(sys);
gen_op Ud1 = Rz(sys, +90)*prop(H, t1dt); // Delay t1 + TPPI propagator
gen_op D = Fm(sys); // Detector to F-
gen_op Upx = lypuls_U(sys, 0., 90.); // Propagator for 90-x pulse
gen_op Uphi[P_cycl], Uphiy; // Propagators for phi pulses
for(int j=0; j<P_cycl; j++)
{
    Uphiy = lypuls_U(sys, 90., P_mix[j]); // Generate phi-y pulse props
    Uphi[j] = Uphiy; // Store this phi-y pulse props
}
```

```
gen_op sigma3, sigman;
```

```
//           Apply Pulse Sequence
```

```
File ecosy; // Declare and open file
ecosy.open("ecosy.dat",
    io_writeonly, a_create);
block_1D t2BLK(t2pts); // Set 1D block for output
gen_op sigma0 = sigma_eq(sys); // Set density matrix equilib.
gen_op sigma1 = evolve(sigma0, Upx); // Apply first (PI/2)x pulse
gen_op sigma2 = sigma1; // Initial sigma2 (t1 = 0)
for(int t1=0; t1<t1pts; t1++) // Loop over all t1 increments
{
    sigma3 = sigman; // Set sigma3 to NULL
    for(int i=0; i<P_cycl; i++) // apply E-COSY phase cycle
        sigma3 += P_det[i]*evolve(sigma2, Uphi[i]); //
    FID(sigma3, D, H, t2dt, t2pts, t2BLK); // acquisition
    Felix(ecosy, t2BLK); // output block: Felix
    evolve_ip(sigma2, Ud1); // evolution next t1 + TPPI
}
ecosy.close(); // Close file
double Om = sys.Omega(); // Spectrometer frequency
Felix2D_params(cout, Om, 2.0*NyqF, // Output Felix parameters
    2*t2pts, offset); // Keep screen nice
cout << "\n\n";
}
```

2.4.3 Discussion

Define Constants: In this program the E-COSY phase cycle is explicitly defined; A 5 step cycle with associated phases and weights are stored in accordance with Figure 0-21.

Read in Parameters: As in prior simulations the spin system is read in from a disk file, the name may be supplied directly when running the program. The user is also asked for the acquisition size, whether to use a strong or weak coupling Hamiltonian, any desired offset, and a Nyquist frequency. Note that because TPPI will be used for phase sensitivity along the t_1/f_1 axis, twice as many points are taken along this axis that along the t_2 axis, and the t_1 incrementation time is half of the dwell time.

Set up Operators, Superoperator: The Hamiltonian is first formulated for either strong or weak coupling depending upon the input choice. Next, the t_1 time increment propagator is formulated. Here, in order to produce an acquisition with TPPI, a 90 degree rotation about the z-axis is added so that each incremented step includes the TPPI phase shift. The detector is set to F- in the next line. A propagators for the 90x pulse is now computed, followed by generation of propagators for all five y-pulses of angle ϕ . Operators are declared, the first for a working density matrix and the second for a null operator.

Apply Pulse Sequence: A file name `ecosal` is created and opened. A 1D data block called `t2BLK` is constructed to temporarily store the FID. Following the pulse sequence diagram, the initial density matrix is set to equilibrium. The next step is to apply a 90 degree pulse (without phase cycling) along the x-axis. This is evolved next by the propagator which accounts for t_1 incrementation as well as TPPI phase adjustments, but it is not done on the first step through the t_1 increments, at $t_1=0$. Rather, this step is done at the end of the loop. The next step is evolution under the final pulse of angle ϕ along the y-axis. The E-COSY phase cycle is implemented at this step as discussed previously. This density matrix, `sigma3`, is then used to simulate an FID which is subsequently sent to the file `ecosal.dat` in Felix format. Again, the last step of the loop is actually the earlier t_1 incrementation. The last program steps close the file, get the spectrometer frequency, output Felix parameters necessary for spectral workup, and clear the screen.

2.4.4 Spin System

For this example we will first treat a 4 spin proton system in order verify that the simulation produces correct results. We shall use a system derived from one used in a previous simulation so that we can readily see how to interpret the E-COSY cross peak pattern.

Initial 3-Spin Proton System for Alternate E-COSY Verification




Figure 0-22 Simple 3 and 4 spin protons systems to be utilized in verification of the alternate E-COSY simulation program. The three spin system is identical to that used in the previous example and the 4 spin system is directly derived from it.

ecosy2.sys

SysName	(2) : ECOSY2	- Name of the Spin System
NSpins	(0) : 4	- Number of Spins in the System
Iso(0)	(2) : 1H	- Spin Isotope Type
Iso(1)	(2) : 1H	- Spin Isotope Type
Iso(2)	(2) : 1H	- Spin Isotope Type
Iso(2)	(2) : 1H	- Spin Isotope Type
PPM(0)	(1) : -.1	- Chemical Shifts in PPM
PPM(1)	(1) : 0.0	- Chemical Shifts in PPM
PPM(2)	(1) : .1	- Chemical Shifts in PPM
J(0,1)	(1) : 8.0	- Coupling Constants in Hz
J(0,2)	(1) : 3.0	- Coupling Constants in Hz
J(1,2)	(1) : 2.5	- Coupling Constants in Hz
Omega	(1) : 500	- 1H Spectrometer Frequency in MHz

2.4.5 Workup

The 2D-data sets produced from this simulation is processed by the identical procedure detailed in the previous example, using an acquisition size of 512 points and the provided macro ecosy1.mac.

2.5 Complimentary E-COSY with Superoperators

2.5.1 Description

As was mentioned at the beginning of this chapter, complimentary E-COSY experiments involve only a change in the weighting factors relative to the cycled pulse phases.

Complimentary E-COSY & E-COSY Pulse Sequence with Phase Cycle




Figure 0-23 The E-COSY sequence used in the first example. The complimentary E-COSY associates the same weights, W , with different pulse phases ϕ . Here ϕ_{CE} is used for complimentary E-COSY and ϕ_E for E-COSY.

We shall not review the mathematical formalism for this simulation as it can be found in the previous example. Rather, we shall make only minor modifications to the GAMMA program implementing E-COSY in the first example to produce a program which will generate either the E-COSY or complimentary E-COSY sequence. The focus will be on understanding the difference between E-COSY and complimentary E-COSY experiments.

2.5.2 Program

```

/* ecosy3.cc *****-c++-
**
**                               Example Program for the GAMMA Library
**
** This program simulates either an E-COSY or complimentary
** E-COSY experiment with complete phase cycling. The pulses
** are taken to be ideal and the effects of relaxation are
** not considered.
**
** See: C. Griesinger, O.W. Sorensen, and R.R. Ernst, JMR,
** 75, 474-492, (1987), "Practical Aspects of the E-COSY
** Technique. Measurement of Scalar Spin-Spin Coupling
** Constants in Peptides". The required phase cycle is
** found on page 477 of this reference.
**
** Note: This program treats only homonuclear spin systems.
**
** Superoperators are used to implement the phase cycle
** which speeds up the computation for small systems
** but this program is cumbersome for large systems.
**
*****

#include "gamma.h"
//                               Define Constants

const P_cycl = 12;                // Phase cycle length
const double P_mix[12] =          // Pulse phase cycle, E-COSY
{0,60,0,60,0,60,0,300,120,300,240,300};
const double PC_mix[12] =         // Pulse phase cycle, CE-COSY
{180,240,180,240,180,240,180,120,300,120,60,120};
//                               Begin Program

main (int argc, char* argv[])
{
    cout << "\n\tSimulation of E-COSY and Complimentary E-COSY\n";
    cout << "\t (TPPI Mode, Hard Pulses, No Relaxation)\n";
//                               Read in Parameters

    String filename;               // Name of spin system file
    query_parameter(argc, argv, 1, // Get filename from command
        "\n\tSpin system filename? ", filename);
    spin_system sys;               // Declare spin system sys
    sys.read(filename);             // Read system from filename
    int t1pts, t2pts;

```

```

    query_parameter(argc, argv, 2, // Get number FID of points
        "\n\tAcquisition Size? ", t2pts);
    t1pts = 2*t2pts;               // Set t1 size for TPPI
    String J;                       // Weak or strong coupling
    query_parameter(argc, argv, 3, // Type of E-COSY desired
        "\n\tWeak or strong coupling (w/s)? ", J);
    String CE;                     // Ask for an offset frequency
    query_parameter(argc, argv, 4, // Choose a Nyquist frequency
        "\n\tComplimentary E-COSY or E-COSY (c/e)? ", CE);
    double offset = query_offset(sys, 0, 1); // Dwell time, quadrature
    double NyqF = query_Nyquist(sys, 0);     // t1 time increment, TPPI
    double t2dt = 1.0/(2.0*NyqF);
    double t1dt = t2dt/2.0;

//                               Set Up Operators, Superoperator

    gen_op H;                       // Set Hamiltonian for
    if(J == "w")                     // Strong or weak coupling
        H = Hcs(sys) + HJw(sys);
    else
        H = Hcs(sys) + HJ(sys);
    gen_op Ud1 = Rz(sys, +90)*prop(H, t1dt); // Delay t1 + TPPI propagator
    gen_op D = Fm(sys);              // Detector to F-
    gen_op Upx = lxpuls_U(sys, 0., 90.); // Propagator for x pulse
    gen_op Upmx = lxpuls_U(sys, 0., -90.); // Propagator for -x pulse
    gen_op U_mix;                   // Temporary mixing propagator
    super_op G_mix;                 // Phase cycle superoperator
    gen_op sigma3;

//                               Construct Mixing, Phase Cycle Superoperator

    double conv = acos(-1)/180;
    double P_det = 0;
    if(CE == "c")                   // This for comp. E-COSY
        cout << "\n\tComplimentary E-COSY Phase Cycle\n";
    else
        cout << "\n\tE-COSY Phase Cycle\n";
    cout << "\n\t# scan               betareference\n";
    for(int i=0; i<P_cycl; i++)
    {
        U_mix = Rz(sys, -P_det);    // Detector phase cycle
        U_mix *= Upmx;              // 3rd 90 Pulse (-x)
        if (CE == "c")              // This for comp. E-COSY
        {
            cout << "\t" << i << "\t" << PC_mix[i]; // Output comp. E-COSY phase
            U_mix = Rz(sys, +PC_mix[i]); // Phase shift pulses 1 & 2
        }
        else                        // This for E-COSY
        {
            cout << "\t" << i << "\t" << P_mix[i]; // Output E-COSY phase
            U_mix = Rz(sys, +P_mix[i]); // Phase shift pulses 1 & 2
        }
    }
    cout << "\t" << P_det << "\n";

```

```
U_mix *= Upx; // 2nd 90 Pulse
U_mix.Op_base(H); // Put in eigenbasis of Ho
G_mix += U_transform(U_mix); // Add to U transform superop
P_det = acos(-cos(P_det*conv))/conv; // Adjust detector phase
}

// Apply Pulse Sequence

File ecosy; // Declare and open file
ecosy.open("ecosy.dat",
io_writeonly, a_create);
block_1D t2BLK(t2pts); // Set 1D block for output
gen_op sigma0 = sigma_eq(sys); // Set density matrix equilib.
gen_op sigma1 = evolve(sigma0, Upx); // Apply first (PI/2)x pulse
gen_op sigma2 = sigma1; // Initial sigma2 (t1 = 0)
for (int t1=0; t1<t1pts; t1++) // Loop over all t1 increments
{
    sigma3 = evolve(sigma2, G_mix); //superop phase cycle
    FID(sigma3,D,H,t2dt,t2pts,t2BLK); //acquisition
    Felix(ecosy, t2BLK); //output block: Felix
    evolve_ip(sigma2, Ud1); //evolution next t1
}
ecosy.close(); // Close file
double Om = sys.Omega(); // Spectrometer frequency
Felix2D_params(cout, Om, 2.0*NyqF,
                2*t2pts, offset); // Output Felix parameters
cout << "\n\n"; // Keep screen nice
}
```

2.5.3 Discussion

Define Constants: Prior to entering the program, the size of the phase cycle (12) and the phase angles for the pulses are set for both the E-COSY and complimentary E-COSY experiment.

Read in Parameters: As in prior simulations the spin system is read in from a disk file, the name may be supplied directly when running the program. The user is also asked for the acquisition size, whether to use a strong or weak coupling Hamiltonian, any desired offset, and a Nyquist frequency. In the program the user is also asked whether E-COSY or complimentary E-COSY is desired. Note that because TPPI will be used for phase sensitivity along the t_1/f_1 axis, twice as many points are taken along this axis than along the t_2 axis, and the t_1 incrementation time is half of the dwell time.


Set up Operators, Superoperator: The Hamiltonian is first formulated for either strong or weak coupling depending upon the input choice. Next, the t_1 time increment propagator is formulated. Here, in order to produce an acquisition with TPPI, a 90 degree rotation about the z-axis is added so that each incremented step includes the TPPI phase shift. The detector is set to F- in the next line. Propagators for the two 90 pulses are now computed, the first along x and the second along -x. A temporary operator is declared for used in summing over the phase cycle. The next line declares a superoperator which will ultimately account for the phase cycle. Finally, a working density matrix is specified.

Construct Mixing, Phase Cycle Superoperator: As outline mathematically in the first example, the entire phase cycle can be removed from the t_1 and t_2 looping by use of a superoperator. In this section, the phase cycle loop is applied and the appropriate superoperator constructed, now the cycle depends upon whether E-COSY or complimentary E-COSY has been specified. The loop goes over the twelve steps of the cycle (defined in the constants section) and at each step the phases are output to the terminal. Initially the detector phase is set to zero. The propagators for each step are then multiplied together in the opposite order of the pulse sequence due to the ordering of the unary *= step. At the end of each cycle, the unitary transform superoperator is summed to account for the mixing propagator. The last step in this section adjusts the detector phase to +/- 180.

Apply Pulse Sequence: Since the superoperator intrinsically contains the phase cycling, the pulse sequence application is quite simple. A file name ecosy is created and opened. A 1D data block called t2BLK is constructed to temporarily store the FID. Following the pulse sequence diagram, the initial density matrix is set to equilibrium. The next step is to apply a 90 degree pulse (without phase cycling) along the x-axis. This is evolved next by the propagator which accounts for t_1 incrementation as well as TPPI phase adjustments, but it is not done on the first step through the t_1 increments, at $t_1=0$. Rather, this step is done at the end of the loop. The next step is evolution under the superoperator which contains all the phase cycling. This density matrix, sigma3, is then used to simulate an FID which is subsequently sent to the file ecosy.dat in Felix format. Again, the last step of the loop is actually the earlier t_1 incrementation. The last program step simply closes the file.

2.5.4 Spin System

For this example we will use the same three spin system as in our first example.



The E-COSY simulation will be run on a number of 3- and 4- spin system which contain protons and coupling constants for this system. These are given below. Only two of the systems are given in full, the others are derived from these two.

ecosy1a.sys

SysName	(2) : ECOSY1A	- Name of the Spin System
NSpins	(0) : 3	- Number of Spins in the System
Iso(0)	(2) : 1H	- Spin Isotope Type
Iso(1)	(2) : 1H	- Spin Isotope Type
Iso(2)	(2) : 1H	- Spin Isotope Type
PPM(0)	(1) : -.1	- Chemical Shifts in PPM
PPM(1)	(1) : 0.0	- Chemical Shifts in PPM
PPM(2)	(1) : .1	- Chemical Shifts in PPM
J(0,1)	(1) : 8.0	- Coupling Constants in Hz
J(0,2)	(1) : 3.0	- Coupling Constants in Hz
J(1,2)	(1) : 2.5	- Coupling Constants in Hz
Omega	(1) : 500	- 1H Spectrometer Frequency in MHz

2.5.5 Workup

The 2D-data sets produced from this simulation is processed by the identical procedure detailed in the previous example, using an acquisition size of 512 points and the provided macro ecosy1.mac.

2.5.6 Results

Simulated Complimentary E-COSY and E-COSY Spectra




Figure 0-24 Comparison of the E-COSY and complimentary E-COSY sequences on a 3-spin proton system. The information content is equivalent, the shifts due to the passive spin(s) are in opposite directions.

2.6 Heteronuclear E-COSY

2.6.1 Description

We shall employ a variant of the pulse sequence used in example 2 of this section to simulation an E-COSY experiment on a heteronuclear spin system. Employment of the sequence shown below is discussed in the literature.¹

Heteronuclear E-COSY Pulse Sequence and Phase Cycle




Figure 0-25 : The homonuclear E-COSY pulse sequence with appropriate phase cycle and weighting factors. For $\psi = 180$, the weights may be replaced by negative values (i.e. subtract FID's).

In this experiment, the angles ϕ and ϕ' have the same phase cycles but they cycle independently. The angles and relative weights are given by the following equations

$$\phi_j = j(\pi/N) + \delta_{N,K}\pi/(2K) \quad j = 0, 1, \dots, 2N-1. \quad (8-5)$$

$$W_j = \frac{(-1)^{j+1}}{\sin^2(\beta_j/2)}. \quad (8-6)$$

the latter used when $N=K$. The values in the figure utilized $N=K=3$. The evolution of the density matrix for the heteronuclear E-COSY sequence is virtually identical with the treatment of the homonuclear E-COSY experiment given in the second example of this Chapter (Alternative E-

1. C. Griesinger, O.W. Sørensen, and R.R. Ernst (1986), *J. Chem. Phys.*, **85**, 6837-6851, "Correlation of connected transitions by two-dimensional NMR spectroscopy". See FIG. 4 on page 6839 as well as equation (43) on page 6847.

COSY pulse sequence). The changes are only in that the pulses are selective for specific isotopes and the last pulse is applied to two different isotopes with two different phases. As such we shall here provide only an outline for the simulation program.

Program Outline for the Homonuclear E-COSY Sequence

Step

Density Matrix

1. Form Equilibrium Density Matrix

$$\sigma_{eq}$$

2. Apply the 1st 90 x-pulse on I

$$\sigma_1$$

3. Evolve for t_1 time

$$\sigma_2(t_1)$$

4. Apply phi pulse on S

$$\sigma_3(t_1)$$

5. Apply phi pulse on I, E-COSY Phase Cycle on I


$$\sigma(t_1, t_2)$$

6. Repeat Steps 4 and 5, E-COSY Phase Cycle on S

7. Acquire FID for this t_1 time

$$\sigma(t_1, t_2)$$

8. Repeat Steps 3-7 for each t_1 time



2.6.2 Program

```
/* ecosy4.cc *****-C++-
**
** Example Program for the GAMMA Library
**
** This program simulates either an E-COSY experiment with
** complete phase cycling. The pulses
** are taken to be ideal
** and the effects of relaxation are not considered.
**
** See: C. Griesinger, O.W. Sorensen, and R.R. Ernst, JCP,
** 85(12), 6837-6852, (1986), "Correlation of Connected
** Transitions by Two-Dimensional NMR Spectroscopy". The required
** required pulse sequence is found on page 6839 of this reference
**
** Note: This program treats heteronuclear spin systems.
**
*****/

#include "gamma.h"

void setup_2D(spin_system& sys, String& S, String& I, int& isoS,
             int& isol, double& offS, double& offI, double& NyqFS,
             double& NyqFI, int& t2pts, int& t1pts=-1)
{
    isoS = query_isotope(sys, S, // Choose an isotope to detect
                        "\n\tIsotope type for detection? ");
    isol = query_isotope(sys, I, // Choose an 2nd channel isotope
                        "\n\tIsotope type for heteroatom? ");
    if(isoS == isol)
    {
        cout << "\n\n\tSorry, The Two Spin Types Are Not"
              << "\n\tAllowed to Be The Same\n\n";
        exit(-1);
    }
    offS = query_offset(sys, isoS); // Ask for an S offset frequency
    NyqFS = query_Nyquist(sys, isoS); // Choose an S Nyquist frequency
    offI = query_offset(sys, isol); // Ask for an I offset frequency
    NyqFI = query_Nyquist(sys, isol); // Choose an I Nyquist frequency
    cout << "\n\tAcquisition Size? "; // Get number FID of points
    cin >> t2pts;
    if(t1pts != -1)
    {
        cout << "\n\tNumber t1 Increments? "; // Get number of t1 points
        cin >> t1pts;
    }
    return;
}
```

```
// Begin Program
main (int argc, char* argv[])
{
    cout << "\n\t\tHeteronuclear E-COSY Simulation\n";
    cout << "\t (TPPI Mode, Hard Pulses, No Relaxation)\n";
    // Read in Parameters
    spin_system sys;
    query_sys(argc, argv, 1, sys); // Get the spin system
    if(sys.homonuclear()) // Only heteronuclear allowed
    {
        cout << "\n\n\tPlease Run One of the Homonuclear"
              << "\t E-COSY Example Programs for This System\n\n";
        exit(-1);
    }
    String Isol, IsoS; // Two isotope types
    int isoS, isol, t1pts, t2pts;
    double offS, offI, NyqFS, NyqFI;
    setup_2D(sys, IsoS, Isol, isoS, isol,
             offS, offI, NyqFS, NyqFI, t2pts);
    double t2dt = 1.0/(2.0*NyqFS); // Dwell time, quadrature
    double t1dt = 1.0/(4.0*NyqFI); // t1 time increment, TPPI
    t1pts = 2*t2pts; // Set t1 size for TPPI
    String J;
    query_parameter(argc, argv, 2, // Weak or strong coupling
                  "\n\tWeak or strong coupling (w/s)? ", J);
    // Set Up The Heteronuclear E-COSY Phase Cycle
    int P_cycl = 6; // Phase cycle length
    double P_mix[P_cycl]; // Pulse phase cycle angles
    double P_det[P_cycl]; // Weighting factors
    double sinBo2, sign = 1.0;
    for(int k=0; k<=P_cycl; k++)
    {
        P_mix[k] = double(k)*60.0 + 30.0;
        sinBo2 = sin(P_mix[k]*PI/360.0); // Angle is phi/2 in radians
        sign *= -1.0;
        P_det[k] = sign/(sinBo2*sinBo2);
    }
    // Set Up Operators, Superoperator
    gen_op H; // Set Hamiltonian for
    if(J == "w") // Strong or weak coupling
        H = Hcs(sys) + HJw(sys);
    else
        H = Hcs(sys) + HJwh(sys);
    gen_op Ud1 = Rz(sys, Isol, +90.0)*prop(H, t1dt); // Delay t1 + TPPI prop
    gen_op D = Fm(sys, IsoS); // Detector to F-
    gen_op Upx = Ixypuls_U(sys, Isol, 0., 90.); // Prop for 90-x pulse
    gen_op UphiS[P_cycl], UphiI[P_cycl], Uphiy; // Prop for phi pulses
```

```

for(int jj=0; jj<P_cycl; jj++)
{
  Uphiy=lxypuls_U(sys,IsoS,90.,P_mix[jjj]); // Phi-y S pulse props
  Uphis[jjj] = Uphiy; // Store pulse prop
  Uphiy=lxypuls_U(sys,IsoI,90.,P_mix[jjj]); // Phi-y I pulse props
  Uphil[jjj] = Uphiy; // Store pulse prop
}
gen_op sigma3, sigman, sigmat;

// Apply Pulse Sequence

File ecosy; // Declare and open file
ecosy.open("ecosy.dat",
io_writeonly, a_create);
block_1D t2BLK(t2pts); // Set 1D block for output
gen_op sigma0 = sigma_eq(sys); // Set density matrix equilib.
gen_op sigma1 = evolve(sigma0, Upx); // Apply first (PI/2)x pulse
gen_op sigma2 = sigma1; // Initial sigma2 (t1 = 0)
for(int t1=0; t1<t1pts; t1++) // Loop over all t1 increments
{
  sigma3 = sigman; // Set sigma3 to NULL
  for(int i=0; i<P_cycl; i++) //apply phase cycle I
  {
    sigmat = P_det[i]*evolve(sigma2,Uphil[i]);
    for(int j=0; j<P_cycl; j++) //apply phase cycle S
      sigma3+=P_det[j]*evolve(sigmat,Uphis[j]);
  }
  FID(sigma3,D,H,t2dt,t2pts,t2BLK); //acquisition
  Felix(ecosy, t2BLK); // output block: Felix
  evolve_ip(sigma2, Ud1); //evolution next t1 + TPPI
}
ecosy.close(); // Close file
double OmS = sys.Omega(IsoS); // Spectrometer frequency
double Oml = sys.Omega(IsoI); // Spectrometer frequency
Felix2D_params(cout,OmS,2.0*NyqFS,2*t2pts, // Output Felix parameters
  offS,Oml,2.0*NyqFI,t1pts,offI,2);
cout << "\n\n"; // Keep screen nice
}

```

2.6.3 Discussion

Initial Function: A short function, `setup_2D`, has been placed at the start of this program which asks for parameters suited for a heteronuclear 2D experiment. The user is here asked for spin types, the acquisition size, any desired offsets, and Nyquist frequencies.

Read in Parameters: The spin system is read in from a disk file and then checked to insure it is indeed heteronuclear. The user is asked for all the appropriate parameters for the heteronuclear experiment in the previously set function, `setup_2D`. Delay times are set from the specified Nyquist frequencies. Note that because TPPI will be used for phase sensitivity along the t_1/f_1 axis, twice as many points are taken along this axis than along the t_2 axis, and the t_1 incrementation time is half of the dwell time. The user is then asked whether to use a strong or weak coupling Hamiltonian.


Set Up The Heteronuclear E-COSY Phase Cycle: Unlike the previous examples in this Chapter, the heteronuclear phase cycle and weighting functions are calculated. This is because the weights are not integer values so they are more accurately obtained by direct computation.

Set up Operators, Superoperator: The Hamiltonian is first formulated for either strong or weak coupling depending upon the input choice. Next, the t_1 time increment propagator is formulated. Here, in order to produce an acquisition with TPPI, a 90 degree rotation about the z-axis is added so that each incremented step includes the TPPI phase shift. The detector is set to F- in the next line. A propagators for the 90x pulse for spins of isotope type I is now computed, followed by generation of propagators for all y-pulses of angle phi. The phi pulses are performed on both I and S channels and may be applied with differing phases so the corresponding propagators are stored separately for each of the isotopes. Operators are declared for a working density matrix, a null opera-

tor, and a scratch operator.

Apply Pulse Sequence: A file name `ecosal` is created and opened. A 1D data block called `t2BLK` is constructed to temporarily store the FID. Following the pulse sequence diagram, the initial density matrix is set to equilibrium. The next step is to apply a 90 degree pulse (without phase cycling) along the x-axis. This is evolved next by the propagator which accounts for t_1 incrementation as well as TPPI phase adjustments, but it is not done on the first step through the t_1 increments, at $t_1=0$. Rather, this step is done at the end of the loop. The next step is evolution under the final pulse(s) of angle ϕ along the y-axis. The E-COSY phase cycle is implemented at this step as discussed previously. Because the ϕ pulse is applied to both channels and phase cycled independently on each a double loop runs through all combinations of phases for I and S. The resulting density matrix, `sigma3`, is then used to simulate an FID which is subsequently sent to the file `ecosal.dat` in Felix format. Again, the last step of the loop is actually the earlier t_1 incrementation. The last program steps close the file, get the spectrometer frequencies, output Felix parameters necessary for spectral workup, and clear the screen.

2.6.4 Spin Systems

Heteronuclear E-COSY Spin Systems for Example Runs**ecosy4a.sys**

SysName (2) : ECOSY4A- Name of the Spin System
NSpins (0) : 3 - Number of Spins in the System
Iso(0) (2) : 1H - Spin Isotope Type
Iso(1) (2) : 1H - Spin Isotope Type
Iso(2) (2) : 19F - Spin Isotope Type
PPM(0) (1) : -.1 - Chemical Shifts in PPM
PPM(1) (1) : 0.0 - Chemical Shifts in PPM
PPM(2) (1) : .1 - Chemical Shifts in PPM
J(0,1) (1) : 20.0 - Coupling Constant JAB in Hz
J(0,2) (1) : 12.0 - Coupling Constant JAC in Hz
J(1,2) (1) : 7.0 - Coupling Constant JBC in Hz
Omega (1) : 500 - 1H Spectrometer Frequency in MHz

ecosy4b.sys

SysName (2) : ECOSY4B- Name of the Spin System
NSpins (0) : 4 - Number of Spins in the System
Iso(0) (2) : 1H - Spin Isotope Type
Iso(1) (2) : 1H - Spin Isotope Type
Iso(2) (2) : 19F - Spin Isotope Type
Iso(3) (2) : 19F - Spin Isotope Type
PPM(0) (1) : -.1 - Chemical Shifts in PPM
PPM(1) (1) : 0.0 - Chemical Shifts in PPM
PPM(2) (1) : -.1 - Chemical Shifts in PPM
PPM(3) (1) : .1 - Chemical Shifts in PPM
J(0,1) (1) : 20.0 - Coupling Constant JAB in Hz
J(0,2) (1) : 12.0 - Coupling Constant JAC in Hz
J(0,3) (1) : 19.0 - Coupling Constant JAD in Hz
J(1,2) (1) : 7.0 - Coupling Constant JBC in Hz
J(1,3) (1) : 5.0 - Coupling Constant JBD in Hz
J(2,3) (1) : 15.0 - Coupling Constant JCD in Hz
Omega (1) : 500 - 1H Spectrometer Frequency in MHz

2.6.5 Results





Figure 0-26 Heteronuclear E-COSY simulations on a 3 and 4 spin system. The coupling constants are displayed on the 3 spin. The solid boxes show the shifts on the cross peaks due to the “passive” proton. The boxes are square with sides of length J_{FH} , where J_{FH} is the scalar coupling to the passive spin.

2.7 E-COSY Phase Cycling


2.7.1 Description

In this chapter we have based our E-COSY simulations on a few basic pulse sequences which are found in the literature.

Equivalent Homonuclear E-COSY Pulse Sequences



Heteronuclear E-COSY Pulse Sequence



The most complicated aspect of implementing the E-COSY simulation is determining which phase cycle is appropriate. The program given here calculates the appropriate phases and weighting factors for the E-COSY experiment.

2.7.2 Program

```

/* ecosypc.cc *****-C+--*
**
** This program calculates the appropriate phase angles and weighting
** functions for an E-COSY experiment. The program runs interactive, asking
** the user to supply the maximum coherence order expected (>2).
**
** See: C. Griesinger, O.W. Sorensen, and R.R. Ernst, JCP, 85(12), 6837-
** 6852, (1986), "Correlation of Connected Transitions by Two-Dimensional
** NMR Spectroscopy".
**
** Note: Although the parameters N, K, B0, and B1 are included in the
** formulas which determine the phase cycle, they are restricted to the
** following values:
**
** N = K = coherence order chosen
** B0 = B1 = 0
**
** If different values are desired, the user will have
** to make minor program modifications and re-compile.
***** */

#include "gamma.h"
main (int argc, char* argv[])

{
// Get the Coherence Order for the E-COSY Phase Cycle

int ORD = 2;
query_parameter(argc, argv, 1,          // Coherence order for cycle
"\n\tCoherence Order (>2): ", ORD);
while(ORD < 3)
{
cout << "\n\tThe coherence order must be > 2";
cout << "\n\tCoherence Order (>2): ";
cin >> ORD;
}
int P_cycl = 2*ORD;          // E-COSY phase cycle length
double B0 = 0.0;
double B1 = 0.0;
int N = ORD;
int K = ORD;

//          Cycle for Homonuclear E-COSY Phase Cycle
//   C. Griesinger, O.W. Sorensen and R.R. Ernst, JMR, 75, 474-492 (1987)
// phi = j * pi/N where j = [0, 2N-1]          [3], page 476
//   j
//       K

//          ---
//          \
// W = 0.5*B + / B cos(p*phi ) - 0.5*B cos(K*phi )
// j   0 --- p   j       K       j
//       p
//       2           2
// B = p /4 + B (p even) ; B = (p -1)/4 + B (p odd)
// p       0           p       1
//       1
// comp
// W   = W
// j   j+N

// Pulse phase cycle angles
// Weighting factors
// Weighting factors, compliment.

double phi_homo[P_cycl], phij;
double W_homo[P_cycl];
double W_Chomo[P_cycl];
double W, Bp, Cpfj;
int eo;
for(int j=0; j<P_cycl; j++)
{
phi_homo[j] = double(j)*180.0/double(N);
phij = phi_homo[j]*PI/180.0;
W = 0.5*B0;
eo = 1;
for(int p=1; p<=K; p++)
{
eo *= -1;
Bp = 0.25*double(p)*double(p);
if(eo > 0)
Bp += B0;
else
Bp += (B1 - 0.25);
Cpfj = cos(double(p)*phij);
W += Bp*Cpfj;
}
if(N==K)
W -= 0.5*Bp*Cpfj;
W_homo[j] = W;
}

// Track p even vs. p odd
// Determine Bp
// Determine cos(p*phi)
// Add to W
// Store weight factor

int jpn;
double phijpn;
for(j=0; j<P_cycl; j++)
{
jpn = j+N;
p hijpn = double(jpn)*180.0/double(N);
phijpn *= (PI/180.0);
W = 0.5*B0;
eo = 1;
for(int p=1; p<=K; p++)
{
[4], page 476

[1a,1b], page 475

[4], page 476

// The phase angles
// Convert phi to radians
// Begin weight factor

// The phase angles
// Convert phi to radians
// Begin weight factor

```



```

eo *= -1; // Track p even vs. p odd
Bp = 0.25*double(p)*double(p); // Determine Bp
if(eo > 0)
    Bp += B0;
else
    Bp += (B1 - 0.25);
Cpfj = cos(double(p)*phijpn); // Determine cos(p*phi)
W += Bp*Cpfj; // Add to W
}
if(N==K)
    W -= 0.5*Bp*Cpfj;
W_Chomo[j] = W; // Store total weight factor
}

// Cycle for Heteronuclear E-COSY Phase Cycle
// C. Griesinger, O.W. Sorensen and R.R. Ernst, JCP, 85, 6837-6852 (1986)
// Equations (47) and (48), both on page 6848

// phi = j * pi/N + del * pi/(2*K) where j = [0, 2N-1]
// j N,K

// K * sin[(K+1)*phi] - (K+1) * sin(K*phi)
// j j K
// W = ----- + del * - sin(K*phi)
// j 2 N,K 2 j
// 4 * sin [phi / 2]
// j

double phi_hetero[P_cycl]; // Pulse phase cycle angles
double W_hetero[P_cycl]; // Weighting factors
double W_Chetero[P_cycl]; // Weighting factors, comp.
for(j=0; j<P_cycl; j++)
{
    phij = double(j)*180.0/double(N); // Calculate the angle phi
    if(N == K) phij += 180.0/double(2*K);
    phi_hetero[j] = phij; // Store phi
    phij *= PI/180.0; // Convert phi to radians
    SKfj = sin(double(K)*phij); // Start weight calculation
    SKp1fj = sin(double(K+1)*phij);
    Sfo2 = sin(phij/2.0);
    W = double(K)*SKp1fj - double(K+1)*SKfj;
    W /= (4.0*Sfo2*Sfo2);
    if(K==N) W += double(K)*SKfj/2.0;
    W_hetero[j] = W; // Store the weighting factor
}

for(j=0; j<P_cycl; j++)
{
    jpn = j+N;

```

```

    phijpn = double(jpn)*180.0/double(N); // Calculate the angle phi
    if(N == K) phijpn += 180.0/double(2*K);
    phijpn *= PI/180.0; // Convert phi to radians
    SKfj = sin(double(K)*phijpn); // Start weight calculation
    SKp1fj = sin(double(K+1)*phijpn);
    Sfo2 = sin(phijpn/2.0);
    W = double(K)*SKp1fj - double(K+1)*SKfj;
    W /= (4.0*Sfo2*Sfo2);
    if(K==N) W += double(K)*SKfj/2.0;
    W_Chetero[j] = W; // Store the weighting factor
}

// Output the E-COSY Phase Cycle
cout << "\n\t\t\t\tE-COSY Phase Cycle\n";
cout << "\t\t\t S.A Smith March 1993\n";
cout << "\t\t\t N = " << N << "\t K = " << K << "\n";
cout << "\n Step\t Homonuclear\t\t\t Heteronuclear";
cout << "\n\t\t phi\t W\t Wcomp"
    << "\t phi\t W\t Wcomp\n";
for(j=0; j<P_cycl; j++)
{
    cout << "\n" << dec(j,10)
        << " " << dtoa(phi_homo[j], 'f', 10, 2)
        << " " << dtoa(W_homo[j], 'f', 10, 4)
        << " " << dtoa(W_Chomo[j], 'f', 10, 4)
        << "\t"
        << " " << dtoa(phi_hetero[j], 'f', 10, 2)
        << " " << dtoa(W_hetero[j], 'f', 10, 4)
        << " " << dtoa(W_Chetero[j], 'f', 10, 4);
}

// Output these values into a FrameMaker Table
matrix mx(P_cycl, 7, 0.0); // Construct a matrix
for(int i=0; i<P_cycl; i++) // Fill the matrix
{
    mx.put(i,i,0);
    mx.put(phi_homo[i],i,1);
    mx.put(W_homo[i],i,2);
    mx.put(W_Chomo[i],i,3);
    mx.put(phi_hetero[i],i,4);
    mx.put(W_hetero[i],i,5);
    mx.put(W_Chetero[i],i,6);
}
FM_Mat_Tbl("ecosy_pc.mif", mx); // Output matrix to FM table
cout << "\n\n"; // Keep screen nice
}

```

2.7.3 Discussion

Get Coherence Order: The program begins by asking the

user which coherence order the E-COSY sequence is to be set for. The minimum allowed value is three. The number of phase cycle steps is determined and the values of N and K are set to the specified order (best sensitivity). The values of B0 and B1 are set to zero, they affect the homonuclear E-COSY cycle.

Cycle for Homonuclear E-COSY: In this section the homonuclear phase cycle and weighting factors are determined. The formulas implemented and the reference are in the program text. The first loop produces the phases and weights whereas the second loop determines the weights for the complimentary E-COSY experiment.

Cycle for Heteronuclear E-COSY: This code is the heteronuclear analog of the previous section.

Output the E-COSY Phase Cycle: The determined cycle is output to the console. The data is for the specified coherence order. Both the hetero- and homo-nuclear cases are output and weights for E-COSY and complementary E-COSY.

Output to FrameMaker: For convenience, the data is placed into a table in FrameMaker compatible format (MIF). First a matrix is filled with the data and then the function FM_Mat_Tbl puts the data into the file `ecosal_pc.mif`. The result of the program can be directly incorporated into any FrameMaker document, such as this one.

2.7.4 Results

Table 3: E-COSY Phase Cycle N=K=3

STEP	ϕ	W	W^{comp}	ϕ	W	W^{comp}
	Homonuclear			HeteroNuclear		
0	0	2	0	30	-3.73205	.267949
1	60	-1.5	0.5	90	0.5	-0.5
2	120	0.5	-1.5	150	-.267949	3.73205
3	180	0	2	210	.267949	-3.73205
4	240	0.5	-1.5	270	-0.5	0.5
5	300	-1.5	0.5	330	3.73205	-.267949

Table 19-1 E-COSY phase cycle suited for systems consisting predominately of 3 spin networks.

Table 20: E-COSY Phase Cycle N=K=4


STEP	ϕ	W	W^{comp}	ϕ	W	W^{comp}
	Homonuclear			HeteroNuclear		
0	0	5	1	22.5	-6.56854	-.259892
1	45	-3.41421	-.585786	67.5	.809957	.361616
2	90	1	1	112.5	-.361616	-.809957
3	135	-.585786	-3.41421	157.5	.259892	6.56854
4	180	1	5	202.5	-.259892	-6.56854
5	225	-.585786	-3.41421	247.5	.361616	.809957
6	270	1	1	292.5	-.809957	-.361616
7	315	-3.41421	-.585786	337.5	6.56854	.259892

Table 19-1 E-COSY phase cycle suited for systems consisting predominately of 4 spin networks, i.e. 4 mutually coupled spins.

2.8 E-COSY with Relaxation


2.8.1 Description

We now demonstrate how one may add in the effects of relaxation when simulating the E-COSY experiment. The pulse sequence and phase cycle will be the one used in our original example.



In the figure, W is a weighting factor, ϕ the phase angle of the first two pulses, and ψ the phase angle of the detector. The evolution of the density matrix for the E-COSY sequence is depicted below, in this instance we shall keep our formalism entirely in Liouville space.

E-COSY Pulse Sequence in Terms of Liouville Space Propagators



The formulae for evolution of the density matrix in accordance with the previous figure are as follows.

$$|\sigma_1(\phi)\rangle = \Gamma_p(x, \phi)\sigma_{eq} = \Gamma_z(\phi)\Gamma_p(x)\Gamma_z^{-1}(\phi)\sigma_{eq} \quad (\text{EQ } 9)$$

$$|\sigma_2(t_1, \phi)\rangle = \Gamma_d(t_1)|\sigma_1(\phi)\rangle \quad (\text{EQ } 10)$$

$$\begin{aligned} |\sigma_3(t_1, \phi)\rangle &= \Gamma_p(-x)\Gamma_p(x, \phi)|\sigma_2(t_1, \phi)\rangle \\ &= \Gamma_p(-x)\Gamma_z(\phi)\Gamma_p(x)\Gamma_z^{-1}(\phi)|\sigma_2(t_1, \phi)\rangle \end{aligned} \quad (\text{EQ 11})$$

$$|\sigma(t_1, t_2, \phi)\rangle = \Gamma_d(t_2)|\sigma_3(t_1, \phi)\rangle \quad (\text{EQ 12})$$

The FID is then computed by performing trace operations with a detection operator, F_- , which is also phase cycled.

$$FID(t_1, t_2, \phi, \psi) = \langle F_-^\dagger(\psi) | \sigma(t_1, t_2, \phi) \rangle \quad (12-1)$$

The last step is that we sum the FID's over the phase cycle

$$FID(t_1, t_2) = \sum_{\phi, \psi}^{\text{cycle}} FID(t_1, t_2, \phi, \psi) \quad (12-2)$$

With GAMMA, we could proceed implementing these equations as outlined and directly perform an E-COSY simulation. Obviously the task would then require some scheme for performing the phase cycle, perhaps involving a complicated algorithm and repetitive computations. Rather than using these first six equations as written, we now shuffle around the formulae and attempt to isolate components of the phase cycle in order to simplify implementation of the E-COSY sequence. From the first two equations we have

$$|\sigma_2(t_1, \phi)\rangle = \Gamma_d(t_1)\Gamma_z(\phi)\Gamma_p(x)\Gamma_z^{-1}(\phi)|\sigma_{eq}\rangle .$$

Rotations about the z-axis cannot affect the equilibrium density matrix so the initial rotation super-operator may be safely removed. Furthermore, rotations about z commute with the time evolution propagators due to $[H_0, F_z] = 0$, so we may switch the order of these operations. The result is

$$|\sigma_2(t_1, \phi)\rangle = \Gamma_z(\phi)\Gamma_d(t_1)\Gamma_p(x)\Gamma_z^{-1}(\phi)|\sigma_{eq}\rangle .$$

We define a new density matrix operators which are independent of the phase angle ϕ ,

$$|\sigma_1\rangle = \Gamma_p(x)|\sigma_{eq}\rangle \quad (12-3)$$

and

$$|\sigma_2(t_1)\rangle = \Gamma_d(t_1)|\sigma_1\rangle , \quad (12-4)$$

so that the previous formula, (EQ 2), becomes simply

$$|\sigma_2(t_1, \phi)\rangle = \Gamma_z(\phi)|\sigma_2(t_1)\rangle . \quad (12-5)$$

A cancellation of the rotation superoperator occurs when the next step is formulated. From (EQ 3),

$$\begin{aligned} |\sigma_3(t_1, \phi)\rangle &= \Gamma_p(-x)\Gamma_z(\phi)\Gamma_p(x)\Gamma_z^{-1}(\phi)|\sigma_2(t_1, \phi)\rangle \\ \sigma_3(t_1, \phi) &= \Gamma(-x)\Gamma_z(\phi)\Gamma_p(x)|\sigma_2(t_1)\rangle \end{aligned} \quad (12-6)$$

Note that the phase cycle is explicitly removed from the t_1 incrementation at this point. We now shift our attention to the last sequence step and the detection, propagation for the t_2 delay followed by multiplication with a detection operator respectively. We have for detection

$$F_-(\psi)\sigma(t_1, t_2, \phi) = R_z(\psi)F_-R_z^{-1}(\psi)\sigma(t_1, t_2, \phi) \quad . \quad (12-7)$$

where we can utilize the relationship $Tr\{AB\} = Tr\{BA\}$,

$$\begin{aligned} \langle F_-^\dagger(\psi)|\sigma(t_1, t_2, \phi)\rangle &= Tr\{R_z(\psi)F_-R_z^{-1}(\psi)\sigma(t_1, t_2, \phi)\} \\ &= Tr\{F_-R_z^{-1}(\psi)\sigma(t_1, t_2, \phi)R_z(\psi)\} \end{aligned}$$

and include the t_2 delay step.

$$\begin{aligned} \langle F_-^\dagger(\psi)|\sigma(t_1, t_2, \phi)\rangle &= \langle F_-^\dagger|\Gamma_z^{-1}(\psi)\Gamma_d(t_2)|\sigma_3(t_1, \phi)\rangle \\ &= \langle F_-^\dagger|\Gamma_d(t_2)\Gamma_z^{-1}(\psi)|\sigma_3(t_1, \phi)\rangle \end{aligned} \quad (12-8)$$

Thus, the phase cycle over ψ has been isolated from the detection and t_2 incrementation. If we again focus on the components which depend upon the phase cycle, we have

$$\Gamma_z^{-1}(\psi)|\sigma_3(t_1, \phi)\rangle = \Gamma_z^{-1}(\psi)\Gamma_p(-x)\Gamma_z(\phi)\Gamma_p(x)|\sigma_2(t_1)\rangle \quad . \quad (12-9)$$

The entire E-COSY pulse sequence phase cycle is contained in the above equation and we can formulate a new propagator which contains it. Letting

$$\Gamma_{mix}(\phi, \psi) = \Gamma_z^{-1}(\psi)\Gamma_p(-x)\Gamma_z(\phi)\Gamma_p(x) \quad , \quad (12-10)$$

The trace equation becomes,

$$\langle F_-^\dagger(\psi)|\sigma(t_1, t_2, \phi)\rangle = \langle F_-^\dagger|\Gamma_d(t_2)\Gamma_{mix}(\phi, \psi)\sigma_2(t_1)\rangle \quad . \quad (12-11)$$

In order to compute our 2-dimensional spectrum we need to sum FID's over each phase cycle according to (0-6)

$$FID(t_1, t_2) = \sum_{\phi, \psi}^{\text{cycle}} FID(t_1, t_2, \phi, \psi)$$

2.8.2 Program

```

/* ecosy5.cc *****-C++-
**
**          Example Program for the GAMMA Library
**
** This program simulates an E-COSY experiment with complete
** phase cycling. The pulses are taken to be ideal and the
** effects of relaxation are included during t1 and t2.
**
** See: C. Griesinger, O.W. Sorensen, and R.R. Ernst, JMR,
** 75, 474-492, (1987), "Practical Aspects of the E.COSY
** Technique. Measurement of Scalar Spin-Spin Coupling
** Constants in Peptides". The required phase cycle is
** found on page 477 of this reference.
**
** Note: This program treats only homonuclear spin systems.
**       Superoperators are used to implement the phase cycle
**       which speeds up the computation for small systems
**       but this program is cumbersome for large systems.
**
***** */

#include "gamma.h"

//          Define Constants

const P_cycl = 12;          // Phase cycle length
const double P_mix[12] =   // Pulse phase cycle

{0,60,0,60,0,60,0,300,120,300,240,300};

//          Begin Program

main (int argc, char* argv[])

{
  cout << "\nSimulation of an E-COSY spectrum in TPPI mode\n";
  cout << "\t (Hard Pulses and Relaxation During Delays)\n";
  //          Read in Spectral Parameters

  sys_dynamic sys;          // Dynamic spin system
  query_system(argc, argv, 1, sys); // Read system from file
  int t1pts, t2pts;
  query_parameter(argc, argv, 2,          // Get number FID of points
    "\nAcquisition Size? ", t2pts);
  t1pts = 2*t2pts;          // Set t1 size for TPPI
  String Isotype = sys.symbol(0); // Get isotope type of 1st spin
  int isoset = query_isotope(sys, Isotype); // Choose an isotope type
  double offset=query_offset(sys,isoset,1); // Ask for an offset frequency
  //          Set Up the Isotropic Hamiltonian

  String J;

```

```

  query_parameter(argc, argv, 3,          // Weak or strong coupling
    "\nWeak or strong coupling (w/s)? ", J);
  gen_op H;          // Set Hamiltonian for
  if (J == "w") H = Hcs(sys) + HJw(sys); // strong or weak coupling
  else H = Hcs(sys) + HJ(sys);
  //          Set Up the Relaxation Superoperator, Liouvillian

  super_op R;
  query_parameter(argc, argv, 4,          // Dipolar relaxation
    "\nInclude Dipolar Relaxation (y/n)? ", J);
  if (J == "y") R += RDD(sys, H);
  query_parameter(argc, argv, 5,          // CSA relaxation
    "\nInclude CSA Relaxation (y/n)? ", J);
  if (J == "y") R += RCC(sys, H);
  query_parameter(argc, argv, 6,          // Dipole-CSA relaxation
    "\nInclude Dipole-CSA Relaxation (y/n)? ", J);
  if (J == "y") R += RDCX(sys, H);
  complex icmplx(0,1);          // z = 0 + 1i
  super_op L = icmplx*Hsuper(H); // L = -i*[Ho, ] (rad/sec)
  L += R;          // Full Liouvillian
  double lwhh = LWWh_DD_max(sys); // Expected dipolar linewidth
  double NyqF =          // Choose a Nyquist frequency
    query_Nyquist(sys,isoset,lwhh);
  double t2dt = 1.0/(2.0*NyqF); // Dwell time, quadrature
  double t1dt = t2dt/2.0; // t1 time increment, TPPI

  //          Set Up Other Necessary Operators, Superoperators

  super_op eLt1 = exp(L, -t1dt); // Delay t1 relaxation superop
  gen_op sigma0 = sigma_eq(sys); // Density matrix at equilib
  set_trace(sigma0, 1.0);
  super_op G_d1 = R_prop(eLt1, sigma0); // Delay t1 relaxation prop
  gen_op RTPPI = Rz(sys,90); // Rotation for TPPI phasing
  RTPPI.Op_base(H); // Put op in Ho eigenbasis
  G_d1 = U_transform(RTPPI)*G_d1; // Delay t1 + TPPI propagator
  gen_op D = Fm(sys); // Detector to F-
  gen_op Upx = lxpuls_U(sys, 0, +90); // Propagator for x pulse
  gen_op Upmx = lxpuls_U(sys, 0, -90); // Propagator for -x pulse
  acquire ACQ(D, L, t2dt); // Prepare for acquisitions
  gen_op U_mix; // Temporary mixing propagator
  super_op G_mix; // Phase cycle superoperator
  gen_op sigma3;

  //          Construct Mixing, Phase Cycle Superoperator

  double conv = acos(-1)/180;
  double P_det = 0;
  cout << "\nConstructing Phase Cycle Superoperator\n";
  cout << "\n\t# scan          betareference\n";
  for (int i=0; i<P_cycl; i++)
  {
    cout << "\t" << i << "\n";
    U_mix = Rz(sys,-P_det); // Detector phase cycle

```

```
U_mix *= Upmx; // 3rd 90 Pulse (-x)
U_mix *= Rz(sys,+P_mix[i]); // Phase shift pulses 1 & 2
U_mix *= Upx; // 2nd 90 Pulse
U_mix.Op_base(H); // Put in eigenbasis of Ho
G_mix += U_transform(U_mix); // Add to U transform supero
P_det = acos(-cos(P_det*conv))/conv; // Adjust detector phase
}
// Apply Pulse Sequence
File ecosy; // Declare and open file
ecosy.open("ecosy.dat",
io_writeonly, a_create);
block_1D t2BLK(t2pts); // Set 1D block for output
gen_op sigma1 = evolve(sigma0, Upx); // Apply first (PI/2)x pulse
gen_op sigma2 = sigma1; // Initial sigma2 (t1 = 0)
for (int t1=0; t1<t1pts; t1++) // Loop over all t1 increments
{
    sigma3 = evolve(sigma2, G_mix); // Superop phase cycle
    ACQ(sigma3,t2BLK); // Collect the FID
    Felix(ecosy, t2BLK); // Output block: Felix
    evolve_ip(sigma2, G_d1); // Evolution next t1
}
ecosy.close(); // Close file
}
```

2.8.3 Discussion

Define Constants: In this program the E-COSY phase cycle are explicitly defined; A 12 step cycle and associated phases in accordance with Figure 19-4.

Read in Parameters: As in prior simulations the spin system is read in from a disk file, the name may be supplied directly when running the program. The user is also asked for the acquisition size, whether to use a strong or weak coupling Hamiltonian, any desired offset, and a Nyquist frequency. Note that because TPPI will be used for phase sensitivity along the t1/f1 axis, twice as many points are taken along this axis that along the t2 axis, and the t1 incrementation time is half of the dwell time.

Set up Operators, Superoperator: The Hamiltonian is first formulated for either strong or weak coupling depending upon the input choice. Next, the t₁ time increment propagator is formulated. Here, in order to produce an acquisition with TPPI, a 90 degree rotation about the z-axis is added so that each incremented step includes the TPPI phase shift. The detector is set to F- in the next line. Propagators for the two 90 pulses are now computed, the first along x and the second along -x. A temporary operator is declared for used in summing over the phase cycle. The next line declares a superoperator which will ultimately account for the phase cycle. Finally, a working density matrix is specified.

Construct Mixing, Phase Cycle Superoperator: As outlined mathematically at the start of this section, the entire phase cycle can be removed from the t1 and t2 looping by use of a superoperator. In this section, the phase cycle loop is applied and the appropriate superoperator constructed. The loop goes over the twelve steps of the cycle (defined in the constants section) and at each step the phases are output to the terminal. Initially the detector phase is set to zero. The propagators for each step are then multi-


plied together in the opposite order of the pulse sequence due to the ordering of the unary *= step. At the end of each cycle, the unitary transform superoperator is summed to account for the mixing propagator. The last step in this section adjusts the detector phase to +/- 180.

Apply Pulse Sequence: Since the superoperator intrinsically contains the phase cycling, the pulse sequence application is quite simple. A file name ecosy is created and opened. A 1D data block called t2BLK is constructed to temporarily store the FID. Following the pulse sequence diagram, the initial density matrix is set to equilibrium. The next step is to apply a 90 degree pulse (without phase cycling) along the x-axis. This is evolved next by the propagator which accounts for t_1 incrementation as well as TPPI phase adjustments, but it is not done on the first step through the t_1 increments, at $t_1=0$. Rather, this step is done at the end of the loop. The next step is evolution under the superoperator which contains all the phase cycling. This density matrix, sigma3, is then used to simulate an FID which is subsequently sent to the file ecosy.dat in Felix format. Again, the last step of the loop is actually the earlier t_1 incrementation. The last program steps close the file, get the spectrometer frequency, output Felix parameters necessary for spectral workup, and clear the screen.

2.8.4 Spin System

For this example we will first treat a the simple 3 spin proton system which was investigated in a previous example.

3-Spin Proton System for Homonuclear E-COSY with Relaxaiton



ecosal5.dsys

SysName (2) : ECOSY1A-	Name of the Spin System
NSpins (0) : 3	- Number of Spins in the System
Iso(0) (2) : 1H	- Spin Isotope Type
Iso(1) (2) : 1H	- Spin Isotope Type
Iso(2) (2) : 1H	- Spin Isotope Type
PPM(0) (1) : -.1	- Chemical Shifts in PPM
PPM(1) (1) : 0.0	- Chemical Shifts in PPM
PPM(2) (1) : .1	- Chemical Shifts in PPM
J(0,1) (1) : 8.0	- Coupling Constants in Hz
J(0,2) (1) : 3.0	- Coupling Constants in Hz
J(1,2) (1) : 2.5	- Coupling Constants in Hz
Omega (1) : 500	- 1H Spectrometer Frequency in MHz

2.8.5 Workup

The 2D-data sets produced from this simulation were processed with the program Felix. Each data set was subjected to the following Felix macro, called `ecosal.mac` (the !comments to the right are not to be included in the macro - these are no longer allowed in the current Felix version.). Note that this macro *assumes that the simulation had an acquisition size of 512 points*.

```
def datfile ecosy          ! simulated data file name "ecosy.dat"
def matfile ecosy          ! matrix file name "ecosy.mat"
def t1max 1024             ! dimension t1 size
def t2max 512              ! dimension t2 size
def total 1024
ty Opening source file "&datfile" .....! comment showing which data file
ty Opening matrix file "&matfile" .....! comment showing which matrix file
cmx                         ! close any existing matrix
cl                           ! close any existing dat files
mat &matfile write          ! open matrix ecosy.mat
ty Filling matrix with data ... ! comment that data processing starting
lb 4                        ! set line broadening for apodization
ph0 90                     ! set zero order phase correction
ph1 0                      ! set first order phase correction
ty Complex-FT of t2 into F2 ... ! comment starting F2 transformations
for row 1 &t1max             ! loop through all t1 rows
re &datfile                 ! read block of ecosy.dat file
si &t2max                   ! set the size
em                          ! exponential multiplication
sb &t2max 90                ! sine bell windowing function
zf &total                   ! zero fill to total size
ft                          ! complex FFT
red                         ! reduce data to real
sto 0 &row                  ! store processed data into matrix
next                        ! loop back for next block
ty Real-FT of t1 into F1 ..... ! comment for t1 workup
for col 1 &total            ! loop through all columns
loa &col 0                  ! retrieve column from matrix
si &t1max                   ! specify the size
em                          ! exponential multiplication
sb &t1max 90                ! sine bell
zf 2048                    ! zero fill (will loose with real transform)
rft                         ! real FFT
ph                           ! apply phase adjustment
red                         ! reduce data to real
sto &col 0                  ! store back into matrix
next                        ! get next column
ty Plotting contours .....   ! comments for processing point
ty zooming cross peaks .....
cmx                         ! clear matrices
mat &matfile read           ! open ecosy.mat for reading only
lim 1 1 1024! set plot limits
lim 2 1 1024
lvl 5e-4                   ! set scaling
cpn 1                      ! contour positive and negative peaks
pen 1                      ! set first pen to pen 1
nl 4                       ! four contour levels
cyc 2                      ! pen cycle
```

```
cp                ! contour plot
ty Processing finished ..... ! comment for processing end
end
```

Prior to running this Felix macro the matrix `ecosity.mat` must be constructed with the following Felix command

```
bld ecosy 2 1024 1024 0      ! real array
```

which construct a 2 dimensional array (1K x 1K) of real numbers. After processing the matrix axes are set with the `rmx` command (used twice, once for each axis).

```
rmx
  1(2)                ! specify the dimension (1 or 2)
  500! spectrometer frequency
  128                  ! spectral width
  3                    ! PPM plotted on axis
  512                  ! reference point
  0                    ! reference point value
  F2(F1)              ! axis label (F2 or F1)
```


The simulated E-COSY spectra were placed into this FrameMaker document in encapsulated Postscript format which resulted from sending the HPGL output of Felix through the filter program `hpgltoeps` provided by FrameMaker. To generate the `hpgl` contour plot file, the Felix commands are (once `cp` produces the screen plot)

```
hdv felix.hpgl      ! hardcopy device is file felix.hpgl
hpm 32               ! hardcopy plot mode is hpgl
hcp                 ! produce the plot
```

2.8.6 Results

Simulated E-COSY Spectra on a 3-Spin Proton System

ecosyl1a.sys




The simulation was repeated for the two 3-spin systems, the input spin systems files are used to label the plots. Each plot was output from Felix, placed into FrameMaker MIF format, rescaled and annotated, and finally placed onto this page. All were performed with weak coupling.

2.9 E-COSY with Phase Cycling

2.9.1 Description

We now move up to a more sophisticated simulation, that of an E- COSY (Exclusive CORrelation SpectroscopY) experiment¹. The pulse sequence and phase cycle are given below.


Sorensen Heteronuclear E-COSY Pulse Sequence



-
1. C. Griesinger, O.W. Sørensen, and R.R. Ernst (1985), *JACS*, **107**, 6394-6396, "Two-Dimensional Correlation of Connected NMR Transitions".
C. Griesinger, O.W. Sørensen, and R.R. Ernst (1986), *J. Chem. Phys.*, **85**, 6837-6851, "Correlation of connected transitions by two-dimensional NMR spectroscopy".
C. Griesinger, O.W. Sørensen, and R.R. Ernst (1987), *JMR*, **75**, 474-492, "Practical Aspects of the E.COSY Technique. Measurement of Scalar Spin-Spin Coupling Constants in Peptides".

In the figure, W is a weighting factor, ϕ the phase angle of the first two pulses, and ψ the phase angle of the detector. As we must employ a phase cycling scheme, we first consider the necessary mathematics for performing such a simulation, ultimately employing superoperators in order to simplify the mathematical formalism and reduce computation time. The evolution of the density matrix for the E-COSY sequence in terms of propagators is depicted below.

E-COSY Sequence in Terms of Hilbert Space Propagators



Although the previous figure does not show the detector phase angle, ψ , it will become explicit when formulating the acquisition. The formulae for evolution of the density matrix in accordance with the previous figure are as follows.

$$\sigma_1(\phi) = U_p(x, \phi) \sigma_{eq} U_p^{-1}(x, \phi) = R_z(\phi) U_p(x) R_z^{-1}(\phi) \sigma_{eq} R_z(\phi) U_p^{-1}(x) R_z^{-1}(\phi) \quad (\text{EQ 13})$$

$$\sigma_2(t_1, \phi) = U_d(t_1) \sigma_1(\phi) U_d^{-1}(t_1) \quad (\text{EQ 14})$$

$$\begin{aligned} \sigma_3(t_1, \phi) &= U_p(-x) U_p(x, \phi) \sigma_2(t_1, \phi) U_p^{-1}(x, \phi) U_p^{-1}(-x) \\ &= U_p(-x) R_z(\phi) U_p(x) R_z^{-1}(\phi) \sigma_2(t_1, \phi) R_z(\phi) U_p^{-1}(x) R_z^{-1}(\phi) U_p^{-1}(-x) \end{aligned} \quad (\text{EQ 15})$$

$$\sigma(t_1, t_2, \phi) = U_d(t_2) \sigma_3(t_1, \phi) U_d^{-1}(t_2) \quad (\text{EQ 16})$$

The FID is then computed by performing trace operations with a detection operator, F_- , which is also phase cycled.

$$FID(t_1, t_2, \phi, \psi) = \text{Tr}\{F_-(\psi) \sigma(t_1, t_2, \phi)\} \quad (\text{EQ 17})$$

The last step is that we sum the FID's over the phase cycle

$$FID(t_1, t_2) = \sum_{\phi, \psi}^{\text{cycle}} FID(t_1, t_2, \phi, \psi) \quad (\text{EQ 18})$$

With GAMMA, we could proceed implementing these equations as outlined and directly perform an E-COSY simulation. Obviously the task would then require some scheme for performing the phase cycle, perhaps involving a complicated algorithm and repetitive computations. Rather than using these first six equations as written, we now shuffle around the formulae and attempt to isolate components of the phase cycle in order to simplify implementation of the E-COSY sequence. From the first two equations we have

$$\sigma_2(t_1, \phi) = U_d(t_1)R_z(\phi)U_p(x)R_z^{-1}(\phi)\sigma_{eq}R_z(\phi)U_p^{-1}(x)R_z^{-1}(\phi)U_d^{-1}(t_1) \quad .$$

Rotations about the z-axis cannot affect the equilibrium density matrix so the inner rotation operators may be safely removed. Furthermore, rotations about z commute with the time evolution propagators due to $[H_0, F_z] = 0$, so we may switch the order of these operations. The result is

$$\sigma_2(t_1, \phi) = R_z(\phi)U_d(t_1)U_p(x)\sigma_{eq}U_p^{-1}(x)U_d^{-1}(t_1)R_z^{-1}(\phi) \quad .$$

We define a new density matrix operators which are independent of the phase angle ϕ ,

$$\sigma_1 = U_p(x)\sigma_{eq}U_p^{-1}(x) \quad (\text{EQ 19})$$

and

$$\sigma_2(t_1) = U_d(t_1)\sigma_1U_d^{-1}(t_1) \quad , \quad (\text{EQ 20})$$

so that the previous formula, (EQ 14), becomes simply

$$\sigma_2(t_1, \phi) = R_z(\phi)\sigma_2(t_1)R_z^{-1}(\phi) \quad . \quad (\text{EQ 21})$$

A cancellation of these rotation operators occurs when the next step is formulated. From (EQ 15),

$$\begin{aligned} \sigma_3(t_1, \phi) &= U_p(-x)R_z(\phi)U_p(x)R_z^{-1}(\phi)\sigma_2(t_1, \phi)R_z(\phi)U_p^{-1}(x)R_z^{-1}(\phi)U_p^{-1}(-x) \\ \sigma_3(t_1, \phi) &= U_p(-x)R_z(\phi)U_p(x)\sigma_2(t_1)U_p^{-1}(x)R_z^{-1}(\phi)U_p^{-1}(-x) \end{aligned} \quad (\text{EQ 22})$$

Note that the phase cycle is explicitly removed from the t_1 incrementation at this point. We now shift our attention to the last sequence step and the detection, propagation for the t_2 delay followed by multiplication with a detection operator respectively. We have for detection

$$F_-(\psi)\sigma(t_1, t_2, \phi) = R_z(\psi)F_-R_z^{-1}(\psi)\sigma(t_1, t_2, \phi) \quad . \quad (\text{EQ 23})$$

where we can utilize the relationship $Tr\{AB\} = Tr\{BA\}$,

$$\begin{aligned} Tr\{F_-(\psi)\sigma(t_1, t_2, \phi)\} &= Tr\{R_z(\psi)F_-R_z^{-1}(\psi)\sigma(t_1, t_2, \phi)\} \\ &= Tr\{F_-R_z^{-1}(\psi)\sigma(t_1, t_2, \phi)R_z(\psi)\} \end{aligned}$$

and include the t_2 delay step.

$$\begin{aligned} Tr\{F_-(\psi)\sigma(t_1, t_2, \phi)\} &= Tr\{F_-R_z^{-1}(\psi)U_d(t_2)\sigma_3(t_1, \phi)U_d^{-1}(t_2)R_z(\psi)\} \\ &= Tr\{F_-U_d(t_2)R_z^{-1}(\psi)\sigma_3(t_1, \phi)R_z(\psi)U_d^{-1}(t_2)\} \end{aligned} \quad . \quad (\text{EQ 24})$$

Thus, the phase cycle over ψ has been isolated from the detection and t_2 incrementation. If we again focus on the components which depend upon the phase cycle, we have

$$\begin{aligned} R_z^{-1}(\psi)\sigma_3(t_1, \phi)R_z(\psi) &= \\ &= R_z^{-1}(\psi)U_p(-x)R_z(\phi)U_p(x)\sigma_2(t_1)U_p^{-1}(x)R_z^{-1}(\phi)U_p^{-1}(-x)R_z(\psi) \end{aligned} \quad . \quad (\text{EQ 25})$$

The entire E-COSY pulse sequence phase cycle is contained in the above equation and we can formulate a new propagator which contains it. Letting

$$U_{mix}(\phi, \psi) = R_z^{-1}(\psi)U_p(-x)R_z(\phi)U_p(x) \quad , \quad (\text{EQ 26})$$

The trace equation becomes,

$$Tr\{F_-(\psi)\sigma(t_1, t_2, \phi)\} = Tr\{F_-U_d(t_2)U_{mix}(\phi, \psi)\sigma_2(t_1)U_{mix}^{-1}(\phi, \psi)U_d^{-1}(t_2)\} \quad . \quad (\text{EQ 27})$$

In order to compute our 2-dimensional spectrum we need to sum FID's over each phase cycle according to (EQ 18)

$$FID(t_1, t_2) = \sum_{\substack{\text{cycle} \\ \phi, \psi}} FID(t_1, t_2, \phi, \psi)$$

and we have a different mixing propagator for each combination of ϕ and ψ . We now employ the power of superoperators in order to simplify this cycle. We form the unitary transformation superoperator equivalent to the mixing propagator. The equation then becomes,


$$Tr\{F_-(\psi)\sigma(t_1, t_2, \phi)\} = Tr\left\{F_-U_d(t_2)[\hat{\Gamma}_{mix}(\phi, \psi)\sigma_2(t_1)]U_d^{-1}(t_2)\right\} \quad . \quad (\text{EQ 28})$$

Now when we sum over the phase cycle we have

$$\begin{aligned}
 FID(t_1, t_2) &= \sum_{\phi, \psi} FID(t_1, t_2, \phi, \psi) \\
 &= \sum_{\phi, \psi}^{cycles} Tr \left\{ F_- U_d(t_2) \left\{ \hat{\Gamma}(\phi, \psi) [U_d(t_1) U_p(x) \sigma_{eq} U_p^{-1}(x) U_d^{-1}(t_1)] \right\} U_d^{-1}(t_2) \right\} \\
 &= Tr \left\{ F_- U_d(t_2) \sum_{\phi, \psi}^{cycles} \hat{\Gamma}(\phi, \psi) [U_d(t_1) U_p(x) \sigma_{eq} U_p^{-1}(x) U_d^{-1}(t_1)] U_d^{-1}(t_2) \right\} \\
 &= \left\{ F_- U_d(t_2) (\hat{\Gamma}_{Cycle} [U_d(t_1) U_p(x) \sigma_{eq} U_p^{-1}(x) U_d^{-1}(t_1)] U_d^{-1}(t_2)) \right\}
 \end{aligned} \tag{EQ 29}$$

The phase cycle needed for the E-COSY simulation is now entirely contained in the superoperator. On a computational level, this means that the phase cycling loop can be removed from the loop over t_1 increments. The code becomes more concise and the computation more efficient¹. The E-COSY pulse sequence in terms of the phase cycle superoperator can be depicted as follows.

E-COSY Sequence Using Superoperator Phase Cycle



The simulation is performed exactly this way. We initially form the equilibrium density matrix, the pulse propagators, the t_1 propagator for t_1 incrementation time, and the superoperator representing the entire E-COSY phase cycle.

1. This is true only for small spin systems, <5 spins. As the system size increases it becomes more difficult for computers to handle superoperators due to the size of the Liouville space.

3 MQF-COSY

3.1 Introduction

The Multiple Quantum Filtered COSY experiment shares the same pulse sequence as the E-COSY sequence but employs different phase cycling. The pulse sequence is shown in the following figure.




Figure 0-27 The MQF-COSY pulse sequence. All three pulses are $P/2$, the phase cycle varies depending upon the coherence order selected.

The required phase cycle depends upon the coherence order to be selected. The acquisition phase ψ typically oscillates between 0 and π , whereas pulse phases increment according to the formula

$$\Delta\phi = \frac{2\pi}{1p} \quad (29-1)$$

where p is the coherence order. A few of these phase cycles are given below

MQF-COSY Phase Cycles

0QF		1QF		2QF		3QF	
ϕ	ψ	ϕ	ψ	ϕ	ψ	ϕ	ψ
0	0	0	0	0	0	0	0
90	0	180	0	90	180	60	180
180	0			180	0	120	0
270	0			270	180	180	180
						240	0
						300	180

Figure 0-28 MQF-COSY pulse sequence phase cycles for coherence orders 0-3.

3.2 MQF-COSY with Superoperators & Felix Output

3.2.1 Description

This example will be a variation on a previous simulation of the E-COSY sequence but employs different phase cycling. The pulse sequence is given below.




Figure 0-29 The MQF-COSY pulse sequence. All three pulses are P/2, the phase cycle varies depending upon the coherence order selected.

The phase cycle needed depends upon the coherence order to be selected. In general, the acquisition phase ψ is oscillating between 0 and π , whereas the pulse phases are incremented according to the formula


where p is the coherence order. A few of these phase cycles are given below.

MQF-COSY Phase Cycles

0QF		1QF		2QF		3QF	
ϕ	ψ	ϕ	ψ	ϕ	ψ	ϕ	ψ
0	0	0	0	0	0	0	0
90	0	180	0	90	180	60	180
180	0			180	0	120	0
270	0			270	180	180	180
						240	0
						300	180

The MQF-COSY simulation will be performed in analogous fashion to the previous E-COSY simulation using a superoperator for the phase cycling. The MQF-COSY pulse sequence in terms of the phase cycle superoperator can be depicted as follows¹.

MQF-COSY Sequence Using Superoperator Phase Cycle



The simulation will be performed exactly this. We initially form the equilibrium density matrix, the pulse propagators, the t_1 propagator for t_1 incrementation time, and the superoperator representing the entire MQF-COSY phase cycle. The program by design will ask for the coherence order to be selected and generate the necessary phase cycle and its corresponding superoperator.

1. For the mathematical basis of this sequence utilizing superoperators see the E-COSY simulation and its description.

2.9.2 Program

```

/* mqfc cosy.cc *****-C++-*****
**
**           Example program for the GAMMA Library
**
** The Program reads a spin system (number of spins, chemical shifts
** and J-coupling constants) from a file and simulates a Multiple Quantum
** Filtered COSY experiment complete with phase cycling.
**
*****/

#include "gamma.h"

//           Define Constants

const int t1pts = 256;           // Number of t1 points
const int t2pts = 128;           // Number of t2 points
const double t1dt = 0.002;       // t1 time increment
const double t2dt = 0.004;       // t2 time increment

//           Begin Program

main (int argc, char* argv[])
{
  cout << "\nSimulation of 2D-MQF-COSY spectra in TPPI mode\n";

  //           Read in Parameters

  String filename;               // Name of spin system file
  query_parameter(argc, argv, 1, // Get filename from command
    "Spin system filename? ", filename);
  spin_system sys;               // Declare spin system sys
  sys.read(filename);            // Read system from filename
  int MQF;
  query_parameter(argc, argv, 2, // Get quantum order desired
    "Quantum order selected? ", MQF);
  MQF = abs(MQF);                // Insure a non-negative value
  int P_cycl;                    // Set up phase cycling
  double P_incr;
  if (MQF)
  {
    P_cycl = 2*MQF;              // MQF > 0
    P_incr = 360/P_cycl;
  }
  else
  {
    P_cycl = 4;                  // Zero quantum filter
    P_incr = 90;                 // Supression up to 2QC
  }
  char J;
  query_parameter(argc, argv, 3, // Weak or strong coupling
    "Weak or strong coupling (w/s)?", J);

  //           Set Up Operators, Superoperator

  gen_op H;                      // Set Hamiltonian for
  if (J == 'w')                   // strong or weak coupling
    H = Hcs(sys) + HJw(sys);
  else
    H = Hcs(sys) + HJ(sys);
  gen_op Ud1 = Rz(sys,+90)*prop(H,t1dt); // Delay t1 + TPPI propagator
  gen_op D = Fm(sys);             // Detector to F-
  gen_op Upx = lxypuls_U(sys, +90, 0); // Propagator for x pulse
  gen_op Upmx = lxypuls_U(sys, -90, 0); // Propagator for -x pulse
  gen_op U_mix;                   // Temporary mixing propagator
  super_op G_mix;                // Phase cycle superoperator
  gen_op sigma3;

  //           Construct Mixing, Phase Cycle Superoperator

  double conv = acos(-1)/180;
  double P_det = 0;
  double sign = -1;
  if(!MQF) sign = 1;              // Constant phase ZQC
  double P_mix = 0;
  cout << "\ntscan #\tbeta\treference\n";
  for ( int i=0; i<P_cycl; i++ )
  {
    cout << "\t" << i << "\t" << P_mix << "\t" << P_det << "\n";
    U_mix = Rz(sys,-P_det);       // Detector phase cycle
    U_mix *= Upmx;                // 3rd 90 Pulse (-x)
    U_mix *= Rz(sys,+P_mix);      // Phase shift pulses 1 & 2
    U_mix *= Upx;                 // 2nd 90 Pulse
    U_mix.Op_base(H);             // Put in eigenbasis of Ho
    G_mix += U_transform(U_mix);  // Add to U transform superop
    P_mix += P_incr;              // Adjust pulse phase
    P_det = acos(sign*cos(P_det*conv))/conv; // Adjust detector phase
  }

  //           Apply Pulse Sequence

  File mqfc cosy;                 // Declare and open file
  mqfc cosy.open("mqfc cosy.dat",
    io_writeonly, a_create);
  block_1D t2BLK(t2pts);         // Set 1D block for output
  gen_op sigma0 = sigma_eq(sys); // Set density matrix equilibrium
  gen_op sigma1 = evolve(sigma0, Upx); // Apply first (PI/2)x pulse
  gen_op sigma2 = sigma1;        // Initial sigma2 (t1 = 0)
  for (int t1=0; t1<t1pts; t1++) // Loop over all t1 increments
  {
    sigma3 = evolve(sigma2, G_mix); // superop phase cycle
    FID(sigma3,D,H,t2dt,t2pts,t2BLK); //acquisition
    Felix(mqfc cosy, t2BLK);       // output block: Felix
    evolve_ip(sigma2, Ud1);        // evolution next t1
  }
  mqfc cosy.close();              // Close file
}

```

3.2.2 Discussion

Define Constants: In this program the dwell times and number of points to use on the t_1 and t_2 axes are defined as constants prior to entering the program. The size of the phase cycle (12) and the phase angles for the pulses are set as well here.

Read in Parameters: As in prior simulations the spin system is read in from a disk file, the name may be supplied directly when running the program. Also queried is whether to use a strong or weak coupling Hamiltonian.

Read in Parameters: As in prior simulations the spin system is read in from a disk file, the name may be supplied directly when running the program. Also queried is whether to use a strong or weak coupling Hamiltonian.

Set up Operators, Superoperator: The Hamiltonian is first formulated for either strong or weak coupling depending upon the input choice. Next, the t_1 time increment propagator is formulated. Here, in order to produce an acquisition with TPPI, a 90 degree rotation about the z-axis is added so that each incremented step includes the TPPI phase shift. The detector is set to F- in the next line. Propagators for the two 90 pulses are now computed, the first along x and the second along -x. A temporary operator is declared for used in summing over the phase cycle. The next line declares a superoperator which will ultimately account for the phase cycle. Finally, a working density matrix is specified.

Construct Mixing, Phase Cycle Superoperator: As outline mathematically at the start of this section, the entire phase cycle can be removed from the t_1 and t_2 looping by use of a superoperator. In this section, the phase cycle loop is applied and the appropriate superoperator constructed. The loop goes over the twelve steps of the cycle (defined in the constants section) and at each step the phases are output to the terminal. Initially the detector phase is set to zero. The propagators for each step are then multiplied together in the opposite order of the pulse sequence due to the ordering of the unary *= step. At the end of each cycle, the unitary transform superoperator is summed to account for the mixing propagator. The last step in this section adjusts the detector phase to +/- 180.

Apply Pulse Sequence: Since the superoperator intrinsically contains the phase cycling, the pulse sequence application is quite simple. A file name ecosy is created and opened. A 1D data block called t2BLK is constructed to temporarily store the FID. Following the pulse sequence diagram, the initial density matrix is set to equilibrium. The next step is to apply a 90 degree pulse (without phase cycling) along the x-axis. This is evolved next by the propagator which accounts for t_1 incrementation as well as TPPI phase adjustments, but it is not done on the first step through the t_1 increments, at $t_1=0$. Rather, this step is done at the end of the loop. The next step is evolution under the superoperator which contains all the phase cycling. This density matrix, sigma3, is then used to simulate an FID which is subsequently sent to the file ecosy.dat in Felix format. Again, the last step of the loop is actually the earlier t_1 incrementation. The last program step simply closes the file.

3.2.3 Example Spin System

For the example we will utilize the a prototypical amino acid proton system pictured below.

4-Spin "Amino Acid" Proton System for MQF-COSY Simulations




Figure 0-30 The MQF-COSY simulation will be run on the above 4 spin system which contain appropriate proton shifts and coupling constants.

The corresponding input spin system (ASCII) file read by the program is presented below.

```
mqfcosy.sys
SysName (2) : system - Name of the Spin System
NSpins (0) : 4 - Number of Spins in the System
Iso(0) (2) : 1H - Spin Isotope Type
Iso(1) (2) : 1H - Spin Isotope Type
Iso(2) (2) : 1H - Spin Isotope Type
Iso(3) (2) : 1H - Spin Isotope Type
PPM(0) (1) : -.420 - Chemical Shifts in PPM
PPM(1) (1) : -.035 - Chemical Shifts in PPM
PPM(2) (1) : .250 - Chemical Shifts in PPM
PPM(3) (1) : .385 - Chemical Shifts in PPM
J(0,1) (1) : 6.8 - Coupling Constants in Hz
J(0,2) (1) : 0.0 - Coupling Constants in Hz
J(0,3) (1) : -.9 - Coupling Constants in Hz
J(1,2) (1) : 9.8 - Coupling Constants in Hz
J(1,3) (1) : 0.0 - Coupling Constants in Hz
J(2,3) (1) : -15 - Coupling Constants in Hz
Omega (1) : 256 - Spectrometer Frequency in MHz (1H based)
```


3.2.4 Workup

The 2D-data sets produced from this simulation were processed with the program Felix. Each data set was subjected to the following Felix macro, called mqfcosy.mac (the !comments to the right are not to be included in the macro - these are no longer allowed in Felix.):

```
def datfile mqfcosy          ! simulated data file name "ecosal.dat"
def matfile mqfcosy          ! matrix file name "ecosal.mat"
def t1max 256                ! dimension t1 size
def t2max 128                ! dimension t2 size
def total 256
ty Opening source file "&datfile" .....! comment showing which data file
ty Opening matrix file "&matfile" .....! comment showing which matrix file
cmx                          ! close any existing matrix
cl                           ! close any existing dat files
mat &matfile write           ! open matrix ecosal.mat
ty Filling matrix with data ...! comment that data processing starting
lb 2                         ! set line broadening for apodization
ph0 90                      ! set zero order phase correction
ph1 0                       ! set first order phase correction
ty Complex-FT of t2 into F2 ...! comment starting F2 transformations
for row 1 &t1max              ! loop through all t1 rows
re &datfile                  ! read block of ecosal.dat file
si &t2max                    ! set the size
em                          ! exponential multiplication
sb &t2max 90                 ! sine bell windowing function
zf &total                    ! zero fill to total size
ft                          ! complex FFT
red                         ! reduce data to real
sto 0 &row                   ! store processed data into matrix
next                        ! loop back for next block
ty Real-FT of t1 into F1 .....! comment for t1 workup
for col 1 &total              ! loop through all columns
loa &col 0                   ! retrieve column from matrix
si &t1max                    ! specify the size
em                          ! exponential multiplication
sb &t1max 90                 ! sine bell
zf 512                      ! zero fill (will loose with real transform)
rft                         ! real FFT
ph                          ! apply phase adjustment
red                         ! reduce data to real
sto &col 0                   ! store back into matrix
next                        ! get next column
ty Plotting contours .....   ! comments for processing point
ty zooming cross peaks .....
cmx                          ! clear matrices
mat &matfile read            ! open ecosal.mat for reading only
lim 1 1 256                  ! set plot limits
lim 2 1 256
lvl 1e-4                     ! set scaling
cpn 1                       ! contour positive and negative peaks
pen 1                       ! set first pen to pen 1
nl 4                        ! four contour levels
cyc 2                       ! pen cycle
cp                          ! contour plot
ty Processing finished ..... ! comment for processing end
```

end

Prior to running this Felix macro the matrix ecosy.mat must be constructed with the following Felix command

 bld ecosy 2 256 256 0 ! real array
which construct a 2 dimensional array (256x256) of real numbers. After processing the matrix axes are set with the rmx command (used twice, once for each axis).


rmx	
1(2)	! specify the dimension (1 or 2)
256	! spectrometer frequency
200	! spectral width
3	! PPM plotted on axis
128	! reference point
0	! reference point value
F2(F1)	! axis label (F2 or F1)

The simulated E-COSY spectra were placed into this FrameMaker document in encapsulated Postscript format which resulted from sending the HPGL output of Felix through the filter program hpgltoeps provided by FrameMaker. To generate the hpgl contour plot file, the Felix commands are (once cp produces the screen plot)

hdv felix.hpgl	! hardcopy device is file felix.hpgl
hpm 32	! hardcopy plot mode is hpgl
hcp	! produce the plot

3.2.5 Results

Simulated MQF-COSY Spectra



4 Mathematical Description

4.1 Overview

All functions provided in the NMR library are described mathematically as well as qualitatively. Since these are generally specific higher level functions we here provide the mathematical description which is at the core of these other descriptions.

4.2 Liouville Equation, Lab. Frame

Because we normally use the density matrix to describe the spin system we now look at its equation of motion.

4.3 Liouville Equation, Interaction Rep.

There are many instances in which the acting Hamiltonian consists of a time independent component as well as a time dependent component. In this case it is usually beneficial to switch into the interaction representation in order to simplify solution of the Liouville equation by effectively removing the time independent Hamiltonian component from the equation. Consider such a situation where the total Hamiltonian is given by

$$\mathbf{H}(t) = \mathbf{H}_I + \mathbf{H}_{II}(t) \quad (29-2)$$

so that the Liouville equation is then

$$i\hbar \frac{d\sigma}{dt} = [\mathbf{H}(t), \sigma] = [\mathbf{H}_I + \mathbf{H}_{II}(t), \sigma] \quad (29-3)$$

Any operator can be switched into the interaction representation by

$$\hat{O} = e^{iH_I t} O e^{-iH_I t} \quad (29-4)$$

where operators in the interaction representation are designated with a hat, $\hat{\cdot}$. The equation of motion is a direct analog of the original Liouville equation

$$i\hbar \frac{d\hat{\sigma}}{dt} = [\hat{\mathbf{H}}_{II}(t), \hat{\sigma}] \quad (29-5)$$

as we will show momentarily. Note that from the switch into the interaction representation we obtain a new equation of motion which does not (explicitly) contain the original time independent component \mathbf{H}_I . There is no guarantee, but hopefully equation (29-5) is easier to solve than (29-3)

and once accomplished the density matrix in the interaction representation switched back into the laboratory frame with inverse of (29-4). Now we show the vality of (29-5) by carrying out the differential. The left hand side of the equation is first expanded.

$$i\frac{d}{dt}\hat{\sigma} = i\frac{d}{dt}\{e^{iH_I t}\sigma e^{-iH_I t}\} = i\left\{(iH_I)e^{iH_I t}\sigma e^{-iH_I t} + e^{iH_I t}\frac{d\sigma}{dt}e^{-iH_I t} - e^{iH_I t}\sigma(iH_I)e^{-iH_I t}\right\}$$

$$i\frac{d}{dt}\hat{\sigma} = ie^{iH_I t}\left\{(iH_I)\sigma + \frac{d\sigma}{dt} - \sigma(iH_I)\right\}e^{-iH_I t} = ie^{iH_I t}\left\{\frac{d\sigma}{dt} + i[H_I, \sigma]\right\}e^{-iH_I t}$$

Equating this to the left hand side of the equation we obtain

$$ie^{iH_I t}\left\{\frac{d\sigma}{dt} + i[H_I, \sigma]\right\}e^{-iH_I t} = [\hat{H}_{II}(t), \hat{\sigma}] \quad .$$

Now we can proceed to isolate the differential of the density matrix in the laboratory frame and expand the double commutator.

$$i\left\{\frac{d\sigma}{dt} + i[H_I, \sigma]\right\} = e^{-iH_I t}[\hat{H}_{II}(t), \hat{\sigma}]e^{iH_I t} = [H_{II}(t), \sigma]$$

$$i\frac{d\sigma}{dt} = [H_{II}(t), \sigma] + [H_I, \sigma] = [H_I + H_{II}(t), \sigma]$$

This being identical with the original Liouville equation, we have show that (29-5) is identical to (29-3).

4.4 Liouville Equation, Rotating Frame

Another common occurance in NMR has the acting Hamiltonian consisting of a time independent component as well as a time dependent component which is smoothly rotating about a specified axis. In such a case it is more beneficial to switch, not into the interaction representation, but into a coordinate system which rotates along with the time dependent Hamiltonian. Such a switch into the rotating frame can simplify solution of the Liouville equation by effectively removing the time dependence from Hamiltonian completely¹. A second use of the rotating frame is NMR to remove some large overall frequency, such as the Larmor frequency, so that all system transitions are seen on a Hertz rather than Megahertz scale². We now consider again the switch into a coordinate system which is rotating about some axis u at a frequency Ω . Generally, any operator can be switched

1. This is exactly what is needed when an external rf-field is applied to the system. Since rf is normally applied in the plane transverse to the direction of the static magnetic field, the xy-plane, it can be seen as oscillating about the z-axis and made to look time independent in the frame rotating about this axis at the rf-field frequency.

into this frame according to

$$\tilde{\sigma} = R_u(\Omega t) \sigma R_u^{-1}(\Omega t) = e^{-iF_u \Omega t} \sigma e^{iF_u \Omega t} \quad (29-6)$$

where an underscore tilde, \sim , is used to designate the rotating frame. What will the Liouville equation look like in this rotating frame? We can determine it directly by expanding the differential of the density matrix with respect to time in the rotating frame.

$$\begin{aligned} i \frac{d}{dt} \tilde{\sigma} &= i \frac{d}{dt} \{ e^{-iF_u \Omega t} \sigma e^{iF_u \Omega t} \} = i \frac{d}{dt} \{ R_u(\Omega t) \sigma R_u^{-1}(\Omega t) \} \\ i \frac{d}{dt} \tilde{\sigma} &= i \left\{ (-iF_u \Omega) e^{-iF_u \Omega t} \sigma e^{iF_u \Omega t} + e^{-iF_u \Omega t} \frac{d\sigma}{dt} e^{iF_u \Omega t} + e^{-iF_u \Omega t} \sigma (iF_u \Omega) e^{iF_u \Omega t} \right\} \\ i \frac{d}{dt} \tilde{\sigma} &= i \left\{ (-iF_u \Omega) \tilde{\sigma} + R_u(\Omega t) \frac{d\sigma}{dt} R_u^{-1}(\Omega t) + \tilde{\sigma} (iF_u \Omega) \right\} \\ i \frac{d}{dt} \tilde{\sigma} &= i \left\{ R_u(\Omega t) \frac{d\sigma}{dt} R_u^{-1}(\Omega t) - i[F_u \Omega, \tilde{\sigma}] \right\} \end{aligned}$$

We know from our original Liouville equation what the differential in the laborator frame is.

$$\begin{aligned} i \frac{d}{dt} \sigma &= R_u(\Omega t) [\mathbf{H}(t), \sigma] R_u^{-1}(\Omega t) + [F_u \Omega, \sigma] \\ i \frac{d}{dt} \tilde{\sigma} &= [\tilde{\mathbf{H}}(t), \tilde{\sigma}] + [F_u \Omega, \tilde{\sigma}] = [\tilde{\mathbf{H}}(t) + F_u \Omega, \tilde{\sigma}] \end{aligned}$$

Thus we obtain the Liouville equation in the rotating frame.

$$i \frac{d}{dt} \tilde{\sigma} = [\tilde{\mathbf{H}}(t) + F_u \Omega, \tilde{\sigma}] \quad (29-7)$$

Like the solution in the interaction representation¹, the hope is that equation (29-7) is easier to solve than (29-3) and once accomplished the density matrix in the rotating frame switched back into the laboratory frame with inverse of (29-6).

2. This mimics using a reference or carrier frequency in an experiment. As a result, all frequencies are seen relative to the carrier. Again, since the Larmor precession of all magnetization occurs about the z-axis (the static field axis) we can switch into a frame rotating about this axis at a frequency near the Larmor frequency to "reference" all measured from that value.

1. There are some authors who call this the interaction representation, I assume because $-F_z \Omega \approx H_0$ in homonuclear systems.

4.5 Liouville Equation, Multiple Rotating Frame

We can take the previous discussion concerning the rotating frame a step farther. In heteronuclear systems it is preferable to switch into a suitable reference frame for each isotope type in the system. That is to say, one wishes to maintain a different reference for each of the differing Larmor frequencies in a system. This presents no complications whatsoever as long as the multiple frames are rotating about the same axis. We simply define selective rotation operator $R_{u, \{i\}}$ which affects only the set of spins $\{i\}$ of a particular isotope type.

$$\tilde{Q}' = R_{u, \{i\}}(\Omega_{\gamma_i} t) O R_{u, \{i\}}^{-1}(\Omega_{\gamma_i} t) = e^{-iF_{u, \{i\}} \Omega_{\gamma_i} t} O e^{iF_{u, \{i\}} \Omega_{\gamma_i} t}$$

In this equation we used the definition

$$F_{u, \{i\}} = \sum_j^{spins} \delta_{\gamma_i \gamma_j} I_{u, j}$$

so that the sum only involves spins sharing the same gyromagnetic ratio. A prime was used to indicate a modified rotating frame and a subscript of γ_i to designate that the frequency Ω_{γ_i} is used for all spins of $\{i\}$ having the same gyromagnetic ratio. The Liouville equation in this selective rotating frame is derived the same as before.

$$i \frac{d}{dt} \tilde{\sigma}' = [\tilde{H}'(t) + F_{u, \{i\}} \Omega_{\gamma_i}, \tilde{\sigma}']$$

There is nothing to prevent us from performing another selective rotation about the axis u at a different frequency and different selectivity. The two applicable equations are

$$i \frac{d}{dt} \tilde{\sigma}'' = [\tilde{H}''(t) + F_{u, \{i\}} \Omega_{\gamma_i} + F_{u, \{j\}} \Omega_{\gamma_j}, \tilde{\sigma}'']$$

and

$$\tilde{Q}'' = R_{u, \{j\}}(\Omega_{\gamma_j} t) \tilde{Q}' R_{u, \{j\}}^{-1}(\Omega_{\gamma_j} t)$$

Since all the selective rotation operators commute (they are applied about the same axis) they have no effect on each other and the order in which they are applied is inconsequential. For as many rotating frames f we care to choose

$$i \frac{d}{dt} \tilde{\sigma} = \left[\tilde{H}(t) + \sum_f F_{u, \{i_f\}} \Omega_{\gamma_{i_f}}, \tilde{\sigma} \right]$$

we can write down the appropriate Liouville equation. The multiple rotating frame is given by

$$\tilde{Q} = \prod_f [R_{u, \{i_f\}}(\Omega_f t)] O \prod_f [R_{u, \{i_f\}}^{-1}(\Omega_f t)] \quad .$$

and the nomenclature gets relatively messy. The end result is that we can work in as many rotating frames as we choose so long as they are rotating about the same axis.

4.6 Liouville Equation Solution, Static Hamiltonian

It is easy to verify that the solution to the Liouville equation under a time independent Hamiltonian is given by

$$\sigma(t) = e^{-iHt} \sigma_o e^{iHt} \quad . \quad (29-8)$$

Here t is the time during which the Hamiltonian has acted on the spin system, σ_o represents the initial state of the system and $\sigma(t)$ the final state. An alternative version of this equation, perhaps more clear, is

$$\sigma(t) = \sigma(t_0 + t_e) = e^{-iHt_e} \sigma(t_0) e^{iHt_e} \quad . \quad (29-9)$$

where the total time $t = t_0 + t_e$, t_0 is some initial time, and t_e is the time over which the density matrix evolves under the effects of the Hamiltonian.

It is common to write this solution in terms of the propagator U ,

$$\sigma(t) = U \sigma_o U^{-1} \quad \text{where} \quad U = e^{-iHt} \quad (29-10)$$

and the time over which the propagator is active is implicit in the propagator itself¹. It is enlightening to examine how the individual elements of $\sigma(t)$ behave under a time independent Hamiltonian. From (29-10) we have,

$$\langle \alpha | \sigma(t_e + t_0) | \beta \rangle = \sum_{\gamma} \sum_{\gamma'} \langle \alpha | U | \gamma \rangle \langle \gamma | \sigma(t_0) | \gamma' \rangle \langle \gamma' | U^{-1} | \beta \rangle \quad .$$

If the density matrix and propagators are in the basis of the static Hamiltonian, H , then the propagator elements will be zero unless they are diagonals.

$$\langle \alpha | \sigma(t_e + t_0) | \beta \rangle = \langle \alpha | U | \alpha \rangle \langle \alpha | \sigma(t_0) | \beta \rangle \langle \beta | U^{-1} | \beta \rangle = \lambda_{\alpha} \lambda_{\beta}^* \langle \alpha | \sigma(t_0) | \beta \rangle$$

1. Under some circumstances it will be very important to know at which time the propagator is applied and the amount of time over which the propagator is active. We will deal with these on a case by case basis.

where λ_α are the eigenvalues of \mathbf{U} which in turn are the exponentials of the eigenvalues of \mathbf{H} .

$$\lambda_\alpha^U = \langle \alpha | \mathbf{U} | \alpha \rangle = \langle \alpha | e^{-i\mathbf{H}t} | \alpha \rangle = e^{-i\omega_\alpha t}$$

In this equation, the frequencies ω_α are associated with the energy levels of \mathbf{H} . The end result is that each element of σ (or any other operator evolving under the effects of \mathbf{H}) oscillates in the complex plane at a specific transition frequency as given by

$$\langle \alpha | \sigma(t_e + t_0) | \beta \rangle = e^{-i\omega_\alpha t} e^{i\omega_\beta t} \langle \alpha | \sigma(t_0) | \beta \rangle = e^{-i\omega_{\alpha\beta} t} \langle \alpha | \sigma(t_0) | \beta \rangle \quad (29-11)$$

where $\omega_{\alpha\beta} = \omega_\alpha - \omega_\beta$.

We can also write down the equivalent to equation (29-10), the solution to the von Neumann equation under a time independent Hamiltonian, in superoperator formalism as¹

$$\sigma(t + t_0) = e^{-i\mathbf{H}t} \sigma(t_0) e^{i\mathbf{H}t} = \mathbf{U} \sigma(t_0) \mathbf{U}^{-1} = \mathbf{U} \sigma(t_0) \quad (29-12)$$

In this context, the superoperator equivalent of the unitary transformation involving the propagator \mathbf{U} , namely \mathbf{U} , is determined from²

$$\mathbf{U} = \mathbf{U} \otimes \mathbf{U}^*$$

Where \mathbf{U}^* is the complex conjugate of the propagator \mathbf{U} and \otimes is a tensor product.

1. See EBW, page 16, equation (2.1.41)
2. See EBW, page 24, equation (2.1.83).

4.7 Example Source Codes

